


Cell Culture Flask

Application

- Tissue culture
- Stem cell research
- Cell biology research
- Cell therapy direction research
- Small-scale cell culture, cell screening
- Biochemistry, medicine, microbiology
- Cytology research in animal and plant and basic medical school
- Genetic engineering, antibody engineering drugs, vaccines, etc. Microplate chromogenic reaction plate use

Main Customers

- College of Life Sciences, College of Animal and Plant, College of Veterinary Medicine, School of Clinical Basic Medicine, College of Agriculture, etc.
- Cell biology laboratory, tissue engineering enterprise, antibody engineering enterprise
- Immunocytotherapy company, stem cell enterprise, medical beauty
- Biological products factory, pharmaceutical factory, dairy product factory, food factory
- Third-party medical testing laboratory, experimental technology service enterprise, etc



T150 U-Shaped Canted Neck Cell Culture Flask



Cell Culture Flask


NEST's Cell Culture Flasks whose cell growth areas ranging from 25cm² to 225cm² are available. These flasks are available as Tissue culture treated or non-treated as well as with a vent cap or plug seal cap to meet your requirements.

● Features

- Made of high clarity, 100% virgin polystyrene.
- Sterilized by E-beam, SAL=10⁻⁶.
- Non-Pyrogenic, DNase/Rnase free.
- Frosted writing and clear graduations.
- Notched bottom for slip free slip stacking.
- Clear lot number for batch traceability.
- Packaged in sterile, zip-sealable bags.

● 3 / 5 - Layer Cell Culture Flask

- Made of high clarity, 100% virgin polystyrene.
- Sterilized by E-beam, SAL=10⁻⁶.
- Non-Pyrogenic, DNase/Rnase free.
- 3-Layer Cell Culture Flask Growth area: 520 cm².
- 5-Layer Cell Culture Flask Growth area: 870 cm².
- Individually packaged in sterile bag.



T150 U-Shaped Canted Neck Cell Culture Flask




Yellow caps for non-treated flask



Vent caps and plug seal caps

Vent caps with 0.22 µm hydrophobic filters to ensure gas exchange without contamination

User guide for seeding cells into a 5-layer Cell Culture Flask



1. Prepare cell suspension of the required concentration, then mix it with medium evenly in a container. A volume of 30-50 mL per layer is recommended.
2. Slowly add the mixed solution into the 5-layer Cell Culture Flask with a serological pipette. To avoid foam or bubbles, it is recommended to set the pipette firmly against the wall, enable the stream to flow along the slope, and save a little liquid in the pipette each time.
Notes: While a 10 ml pipette can disperse the medium at the bottom, a 25 ml pipette will only reach up to the NEST mark to disperse the medium.
3. Position the Multi-layer Flask upright with the NEST mark facing you, tilt it 45° clockwise and let stand in this position for a while to level the liquid in each layer.
4. Gently lay it flat onto the workbench with NEST mark facing upwards.
5. Gently shake it from side to side to distribute cells evenly onto culture surfaces.
Notes: be careful to shake gently to avoid foam or bubbles and spilling liquid from each layer.
6. Transfer the flask to the incubator for incubation

● Culture medium removal


1. Aspiration
Tilt the flask 45° clockwise with the NEST mark facing you. count-er-clock wise to a 45° angle while inverting the Multi-Flask toward you. Then reach the serological pipette into the bottom for fully aspiration.
2. Pouring
Tilt the flask 45° counterclockwise with the NEST mark facing you, pour the spent media from the flask.
Tips: A NEST 10mL serological pipette is suggested for fully aspiration

● Cell harvesting

1. Rinse off the residual serum with buffer, add digestion solution (≥ 5 mL per layer) and mix evenly. Then, follow Steps 3-4 to distribute to dissociating reagent to each layer.
2. Let stand for 2 min, then neutralize and mix with inactivating solution following steps 3-4. Gently swirl to dislodge cells completely.
3. Transfer the solution in a centrifuge tube or other containers by aspiration or pouring.
4. Rinse the flask with buffer for three times, then transfer the buffer into the centrifuge tube for passage and counting.
Tips: Search "NEST Multi-layer Flask" video on YouTube(@nestwuxi4075).

● Cell Culture Flask

| Cell Growth Area (cm ²) | 最大刻度容量(mL) | Rcommended Medium Volume (mL) | Size(mm) | | | Packing | | TC Treated | | Non-Treated | |
|-------------------------------------|------------|-------------------------------|----------|--------|--------|---------|-------|----------------|-----------|----------------|-----------|
| | | | 底高 | Width | Length | /Pack | /Case | Plug seal caps | Vent caps | Plug seal caps | Vent caps |
| 25 | 30 | 5-7.5 | 25.9 | 53.8 | 97.13 | 10 | 20 | 707001 | 707003 | 707011 | 707013 |
| 75 | 225 | 15-22.5 | 35.7 | 89.56 | 160.01 | 5 | 20 | 708001 | 708003 | 708011 | 708013 |
| 150 | 375 | 30-45 | 40.3 | 110.75 | 203 | 5 | 8 | 720001 | 720003 | 720011 | 720013 |
| 175 | 400 | 35-52.5 | 39.1 | 120.51 | 217.9 | 5 | 8 | 709001 | 709003 | 709011 | 709013 |
| 225 | 700 | 45-67.5 | 46.05 | 137 | 238.5 | 5 | 5 | 721001 | 721003 | 721011 | 721013 |
| 520 (3层瓶) | 50 | 60-100 | 60.1 | 120.5 | 203 | 1 | 12 | 731301 | 731302 | / | / |
| 870 (5层瓶) | 50 | 100-150 | 84.3 | 120.5 | 203.6 | 1 | 8 | 731001 | 731002 | / | / |



Precautions:

It is important to handle multi-layer Cell Culture Flasks with caution to avoid the formation of bubbles. The presence of bubbles can lead to the creation of siphon bridge at the baffle, resulting in the upper layer of culture medium flowing down to the bottom.