HandyStep® electronic

Mode d'emploi

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Operating Manual

Please read carefully before use, and follow all operating and safety instructions!

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The device named below fulfills the relevant fundamental requirements of the EC directives and standards listed. In the case of unauthorized modifications to the device, the declaration becomes invalid.

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Safety instructions

This instrument may sometimes be used with hazardous materials, operations, and equipment. To address all of the safety problems associated with its use in such applications is beyond the scope of this manual. It is the responsibility of the user of this instrument to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

1. Observe general instructions for hazard prevention and safety instructions; e.g., wear protective clothing, eye protection and gloves.
2. Observe all safety precautions provided by reagent manufacturers.
If in doubt, contact the manufacturer or supplier.
3. Every user must be acquainted with the operating manual before operation.
4. Use the instrument only for dispensing liquids, with strict regard to the defined Operating Exclusions and Limitations (see page 7, 8).
5. Never use force on the instrument.

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Safety instructions

- continued -

6. Operate only with tip attached (see page 16).
7. Avoid touching the tip orifice when working with hazardous samples.
8. Always work in a way which neither endangers the user nor any other person. Dispense samples into a suitable container, avoiding splashing.
9. Only press the STEP key after ensuring that dispensing the sample will not result in any danger. Dispense into a suitable vessel.
10. **When a tip is inserted and locked the plunger can move inside the tip. If the tip was previously used, make sure that no sample was left inside, or that the tip is directed to a receiving vessel.**
11. Do not use the instrument in potentially explosive environments or with explosive samples.

Note:

Attention: Risk of injury!
Do not touch the plunger seat inside the instrument!

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Safety instructions

- continued -

12. Only use the AC adapter indoors. Protect from humidity.
13. To charge the NiMH battery pack use only the original AC adapter in combination with the original charging dock.
14. Do not replace the original battery pack with other battery packs or cells. Use only original manufacturer's replacement battery packs.
15. Only dispose of completely discharged battery packs. Observe relevant disposal regulations for your community.
16. Use only original manufacturer's accessories and spare parts.
17. In case of malfunction, immediately stop dispensing. Consult the "Troubleshooting" section and contact the manufacturer if needed.

Warning!

Improper use of the instrument or the batteries (short circuit, mechanical damage, overheating, incorrect AC adapter etc.) may, in extreme cases, lead to explosion of the batteries.

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Purpose and limitations

The HandyStep® electronic is a microprocessor controlled battery-operated repetitive pipette. The instrument automatically recognizes new PLASTIBRAND® PD-Tips, which are encoded for tip volume. This allows fastest possible volume selection. If you want to use compatible tips from other manufacturers, you must set the volume manually.

From now on, the term "tips" refers to PD-Tips with size encoding and compatible tips.

The instrument is designed for quick and simple repetitive dispensing of samples, observing the following limitations:

- +15 °C to +40 °C (59 °F to 104 °F)
(instrument and reagents – other temperatures on request)
- vapor pressure up to 500 mbar
- viscosity:
 - 20 mPas with 50 ml PD-Tips
 - 260 mPas with 5 ml PD-Tips
 - 977 mPas with 1.25 ml PD-Tips

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Operating exclusions

During proper operation, the dispensed sample only comes in contact with the tip and not with the HandyStep® electronic.

Never dispense samples which attack polypropylene, polyethylene or polycarbonate (housing). Avoid aggressive fumes (danger of corrosion)!

Do not use aggressive cleaning agents to clean the instrument.

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Note:

The instrument and the tips are not autoclavable (121 °C/250 °F). See sterile tips on page 44.

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

Initial operation

Insert battery pack

1. Open the battery compartment cover and insert battery pack, ensuring that the electrical contacts on the battery pack connect with those of the HandyStep® electronic.
2. Replace the battery compartment cover.

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
Charge battery pack

The battery pack can be charged either inside the HandyStep® electronic or in the charging dock.

Before first use charge the battery pack inside the HandyStep® electronic in the original charging dock for at least 3.5 hours. After 2 - 3 further charging cycles inside the instrument, the full capacity is reached.

Charge indicator is most accurate when battery pack is charged 3.5 hours in the instrument.

Removal of battery pack resets indicator. Upon reinsertion battery-charge frame flashes to indicate that charge status is approximate.

Note:

Stored data and settings are retained when the battery pack is completely discharged or during replacement.

Remove the battery pack when not in regular or daily use.

One single tone indicates that the battery pack has to be charged or replaced. The battery symbol will flash on the display permanently.

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The charging time for a nearly empty battery pack mounted in the instrument is approximately 3.5 hours.

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Charging inside the HandyStep® electronic

Attention:

Before charging, make sure that the AC adapter is appropriate for your line voltage.

1. Plug the AC adapter connector into the side of the charging dock.
2. Place the instrument with the battery pack into the charging dock. A tone indicates proper positioning of the instrument.
3. A green light on the charging dock indicates that the instrument is charging.

Attention:

Upon each battery pack insertion, charge-indicator frame flashes until battery pack is charged for at least 3.5 hours inside instrument.

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Charging the battery pack in the charging dock

The charging dock allows you to charge a separate battery pack while the HandyStep® electronic is in use.

1. Plug the AC adapter connector into the side of the charging dock.
2. Insert the battery pack into the central opening of the charging dock, making sure the battery pack contacts connect with the charging dock contacts.
3. A green light on the charging dock indicates that the battery pack is charging.

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The charging time for a nearly empty battery pack in the charging dock is approximately 2.5 hours.

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Replace battery pack

A fully charged battery pack provides enough power for approximately 2.5 hours of continuous use with a 12.5 ml PD-tip filled with a sample with similar viscosity and density to water.

- A tone indicates the need for recharging.
- The battery symbol will flash on the display.

This indicates the need to recharge the HandyStep® electronic or to insert a spare battery pack.

Before removing the battery, make sure the tip is completely empty.

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Switching on/off

- The HandyStep® electronic can be switched on by inserting and locking a tip, or by pressing the ENTER key.
- When the locking lever is in the open position, the display shows "CLO" (Close).
- Display: **rEF** (reference run): Press the STEP-key to start reference run.
rEF is displayed after every replacement of the battery pack.

Auto Power Off

- The HandyStep® electronic is factory-set to switch off after 3 minutes of inactivity.
- To switch off manually, press the ENTER key briefly

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
Inserting tips

Attention:

Risk of injury! Do not touch the plunger seat inside the instrument!

1. Push the locking lever to its upper stop.
2. Slide a tip laterally into the HandyStep® electronic.
3. Push the locking lever down until you hear a distinct clicking (2 steps).
4. If using encoded PD-Tips from BRAND, the PD-Tip capacity and the word "**TIP**" appear on the display for 3 sec.
5. Then the partial volume to be dispensed in each step is displayed. The default setting is 10 % of the maximum tip capacity.

If using compatible tips (not PD-Tips) you must enter the capacity of the tip yourself (see page 31).

Note:

When using 25 ml and 50 ml tips, the supplied adapter must be used. It can be connected to the tips using a twist lock and can be removed after use.

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
- If the locking lever is not fully engaged, horizontal bars are shown on the display.
- If the locking lever was fully closed but no tip was inserted, "**no TIP**" appears on the display.

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
Removing tips

1. Hold the HandyStep® electronic over a suitable container.
2. Completely discharge the tip by continuously pressing the STEP key (releasing the residual stroke lock), or by using the "E" mode (see page 23).
3. Push the locking lever to its upper stop.
4. Remove the tip.

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Note:

After dispensing highly viscous or crystallizing sample, the tip may wear and require replacement

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Functions of the HandyStep® electronic

The instrument performs the following functions:

1. Dispensing **DISP** (see page 22): Default mode after switching on. A user-defined volume can be dispensed repeatedly, depending on the capacity of the tip and the filling volume.

2. Automatic Dispensing **AUTO-DISP** (see page 24): Identical to DISP-mode, but with automatic dispensing by continuously pressing the STEP key. Releasing the STEP key interrupts the process.

In the AUTO-DISP-mode the instrument has a "**Learn**" function: After pipetting at least three steps, the HandyStep® electronic calculates the average dispensing interval, and will repeatedly dispense at that interval as long as the STEP key is pressed.

3. Pipetting **PIP** (see page 26): In PIP-mode, the HandyStep® electronic functions as a positive displacement non-repeating pipette, particularly useful for volatile, foamy or high viscous materials. A preselected volume is aspirated and dispensed (Positive displacement principle).

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Aspirating sample

Note:

Default setting for the HandyStep® electronic is DISP mode (see page 22).

- The arrow on the display points upward Δ ("Aspirate").
- Push the STEP key once. The tip will be filled completely.

Note:

To interrupt the filling process, press the STEP key. To continue filling, press the Δ key first and hold down the STEP key in succession. The step count on the display will increase.

- During filling, the number of remaining steps is continuously displayed.
The arrow on the display points downward ∇ ("Dispense") when the aspirating is stopped.

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Attention:

After aspirating liquid, the HandyStep® electronic will automatically dispense a small amount of sample in order to compensate for gear backlash.

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DISP-mode

- The arrow on the display points downward ∇ ("Dispense").
- Press the STEP key to dispense the selected step volume (see "Adjusting STEP volume" page 27).
- The display shows the remaining number of steps.
- After dispensing the last full step, a small amount of sample will be left in the tip. To discharge this volume, hold down the STEP key.
- Remove tip (see page 18).

Note:

The tip must be completely empty before you remove it.

Attention:

Dispensing can be interrupted at any time by pressing the Δ key.

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

⇒ siehe Seite 26.

Note:

To empty the tip at any time, press the ∇ key and then hold down the STEP key. The letter E appears on the display and the arrow pointing downward flashes. This process can be interrupted at any time by pressing the Δ key.

A partially filled tip can be completely filled by pressing and holding, in sequence, the Δ key and the STEP key.

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AUTO-DISP-mode

1. To activate the AUTO-DISP-mode, briefly press the MENU key twice.

Note:

The AUTO-DISP-mode must be actively selected.

2. Select the AUTO-DISP-mode with the $\nabla\Delta$ keys (ON/OFF).
3. Confirm by pressing ENTER.
4. "Learn" feature: After pipetting at least three steps, the HandyStep® electronic automatically determines the average time interval between steps (page 19).
5. Hold down the STEP key to dispense the selected volume (page 27) at the calculated intervals.

Note:

Intervals can range from 0.1 to 5 seconds.

6. To interrupt dispensing, release the STEP key. The calculated step interval remains in memory.
7. In AUTO-DISP-mode, the "Learn" process can be restarted at any time.

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PIP-mode

(Reverse pipetting,
positive displacement principle)

1. Select the PIP-mode by holding down the MENU key (> 3s), and press it briefly (page 30).

Make sure that plunger is at lowest position!

2. Select a volume (page 27).

3. To aspirate the selected volume, press the STEP key. A volume slightly larger than the selected volume will be drawn into the tip.

4. To dispense the volume, press the STEP key again. Some residual volume will remain in the tip.

5. To aspirate a new sample, press the STEP key again.

6. When you have completed pipetting, discharge the residual stroke by holding down the STEP key.

Attention!

To interrupt the pipetting process at any time, press the Δ key.

Note:

To change from PIP-mode to AUTO-DISP-mode, you have to switch to DISP-mode first.

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
Adjusting STEP volume

Note:

After the tip size has been recognized, the dispensing volume is shown on the display. The default step volume is 10 % of the nominal capacity.

1. Press the MENU key once. The volume (VOL) is shown on the display.
2. Use the $\nabla\Delta$ keys to select the appropriate volume. Holding down the $\nabla\Delta$ keys will accelerate the rate of change.
3. Confirm the volume selection by pressing ENTER.

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

Adjusting aspiration and dispensing speed

1. Hold down the MENU key (> 3s) to adjust the aspiration speed. The arrow for selection of aspiration speed lights up on the display. Choose between three different aspiration speeds.

Note:

Lowest speed:
Display shows outline of arrow with one bar.
Highest speed:
Display shows outline of arrow with three bars.

2. Use the $\nabla\Delta$ keys to select aspiration speed.
3. Confirm aspiration speed setting by pressing ENTER.
4. The down arrow appears on the display.
5. Use the $\nabla\Delta$ keys to select dispensing speed.
6. Confirm settings by pressing ENTER. You are returned to the main menu.

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Use of compatible tips

Compatible dispenser tips which are not manufactured by BRAND are not encoded for size. These tips can still be used with the HandyStep® electronic ("Insert tips" page 16). After insertion of one of these tips, the display will automatically go to tip selection mode.

By default, the 5 ml tip is displayed flashing regardless of the inserted tip. 12 possible tip sizes are available.

1. Use the $\nabla\Delta$ keys to select appropriate tip size.
2. Confirm setting by pressing ENTER (page 27).

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Residual stroke

A lock prevents the residual stroke from being dispensed as a regular dispensed or pipetted volume. The symbol Δ indicates this remaining volume. It can be discharged by holding the STEP key.

Non-contact dispensing

When tips with volumes of 5 ml or larger are selected, the default dispensing speed adjusts to allow accurate dispensing without contacting the walls of vessel (free-jet dispensing). This can be overridden by selecting a different speed (page 28).

Maintenance and cleaning

The HandyStep® electronic is calibrated at the factory and is maintenance-free. If required, clean the instrument with a cloth dampened with water or a mild soap solution. The HandyStep® electronic must not be dismantled.

Default settings

To revert the HandyStep® electronic to its factory default settings, simultaneously press the MENU and ENTER keys for > 3 seconds, with the plunger in the lowest position.

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Troubleshooting

Problem	Error message displayed	Possible cause	Corrective action
Tip dripping	No message	Leaking tip	Replace tip
Faint display; reduced speed	Battery symbol flashing	Battery low	Recharge battery pack or replace with charged battery pack.
Instrument not responding	ERR 001	<ul style="list-style-type: none">- Locking lever not engaged during aspirating or dispensing- Locking lever open at lower third of the tip	<ol style="list-style-type: none">1. Press ENTER to confirm error. Follow MENU-instructions to avoid instrument damage! Display: CLO (Close)2. Engage the locking lever Display: rEF (Reference run)3. Press the STEP key <p>Attention: Liquid remaining in tip will be discharged during the reference run (rEF). Discharge sample into a suitable vessel (see page 5, point 10).</p>
Instrument not responding	ERR 002	<ul style="list-style-type: none">Motor blocked- Battery low- Viscosity is too high (s. page 7)- Tip orifice obstructed	<ol style="list-style-type: none">1. Press ENTER to confirm error. Follow MENU-instructions to avoid instrument damage! - Recharge or replace battery pack - Reduce aspirating/dispensing speed Display: rEF (Reference run)2. Press the STEP key <p>Attention: Liquid remaining in tip will be discharged during the reference run (rEF). Discharge sample into a suitable vessel (see pages 5, point 10).</p>
Instrument not responding	ERR TIP	Locking lever opened too early within the lower third of the tip	<ol style="list-style-type: none">1. Press ENTER to confirm error. Follow MENU-instructions to avoid instrument damage!2. Display: OPE (Open) Push the locking lever to its upper stop.3. Display alternating between out TIP and CLO. Remove tip and close locking lever (without tip).

			<ol style="list-style-type: none"> 4. HandyStep® electronic automatically carries out a reference run. 5. Display: no TIP Insert new tip (see page 16).
Instrument not responding	ERR 003	Electronic fault or tip connector not at lowest position when inserting new tip	<ol style="list-style-type: none"> 1. Press ENTER to confirm error. 2. Follow MENU-instructions: <ul style="list-style-type: none"> - Display rEF (reference run) - Press STEP-Key <p>Attention: Liquid remaining in tip will be discharged during reference run (rEF). Discharge sample into suitable vessel (see page 5, point 10).</p> <p>If error message is still shown on the display remove battery pack and insert it again (reset). If this does not solve the problem send instrument back for repair.</p>
No tone, no light when placing the HandyStep® electronic into the charging dock	ACC	<ul style="list-style-type: none"> - Instrument contains no battery - Battery pack inserted improperly - Battery pack faulty 	<ul style="list-style-type: none"> - Insert battery pack. - Make sure battery contacts are aligned with contacts in the instrument. - Replace battery pack.
Charge-indicator is flashing	No message	Removal or insertion of a new battery pack	Charge battery pack for at least 3.5 hours inside the HandyStep® electronic.
Dispensing process interrupted. Attention! Volume not completely discharged	ERR VOL	<p>Dispensing process interrupted by</p> <ul style="list-style-type: none"> - Pressing the Δ key during dispensing - electronic faults 	<ol style="list-style-type: none"> 1. Press ENTER to confirm error. 2. Restart dispensing by pressing the STEP key.

Note:

If other error messages are displayed, remove battery pack and insert it again.

If this does not solve the problem, send instrument in for repair.

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Technical data

Battery pack:

Rechargeable Nickel Metal Hydride
battery pack:
NiMH, 4 x AAA cells,
4.8 V/700 mAh

Charging times for empty battery pack:

in battery compartment of charging
dock: approx. 2.5 h

inside the HandyStep® electronic:
approx. 3.5 h

When battery capacity starts to
diminish, recharge battery pack in the
battery compartment of the charging
dock.

Battery pack capacity:

A full charge gives approx. 2.5 hours
of continuous work with a 12.5 ml
PD-Tip.

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Return for repair

Attention:

Transporting of hazardous materials without a permit is a violation of federal law.

BrandTech Scientific, Inc. will not accept instruments that are not appropriately cleaned and decontaminated.

→ Therefore, contact BrandTech Scientific, Inc. and obtain return authorization **before** sending your instrument for service.

Note:

No warranty claims can be considered if repairs other than those described have been made or attempted.

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Warranty

As provided by law, any and all warranties are null and void if the product has been misused, modified or repaired by unauthorized personnel, if the defects are caused by negligence (non-observance of operating manual, lack of maintenance) or by normal wear and tear. Use only original manufacturer's accessory/spare parts.

For conditions and extent of warranty refer to our warranty card enclosed with the product.

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Subject to technical modification without notice. Errors excepted.

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BrandTech Scientific, Inc. no acepta aparatos que no hayan sido debidamente limpiados y descontaminados.