

SPECIMEN ADAPTOR


SPECIMEN ADAPTOR

	EN	BG	HR	CZ	DA	NL
1	SOCK	РЪКАВ	ŽLIJEB ZA ISTAKANJE	TAMPON	SOK	SOK
2	PATIENT PORT	ПОРТ НА ПАЦИЕНТА	PRIKLUČAK ZA PACIJENTA	PORT PACIENTA	PATIENTPORT	PATIËNTPOORT
3	CAP	ЗАПУШАЛКА	POKLOPAC	VÍČKO	HÆTTE	DOP

HI-FLOW COVER

HYDROPHOBIC COVER

	EN	BG	HR	CZ	DA	NL
1	PATIENT	ПАЦИЕНТ	PACIJENT	PACIENT	PATIENT	PATIËNT
2	VACUUM	ВАКУУМ	VAKUUM	PODTLAK	VAKUUM	VACUÛM
3	TANDEM	ДВОЕН	TANDEM	TANDEM	TANDEM	TANDEM
4	POUR SPOUT	ДОЗАТОР	NASTAVAK ZA IZLIJEVANJE	VYLÉVACÍ OTVOR	TAPSNUDE	SCHENKTUIT

SPECIMEN ADAPTOR

INTENDED USE: *Intended for connection to suction canisters to allow for the collection of specimens removed during surgical or medical procedures.*

DIRECTIONS FOR USE - ASSEMBLY

Note: *Begin directions only after properly assembling suction canister.*

1. Insert Sock into Pour Spout.
2. Inspect that Specimen Adaptor is sealed to Canister Cover.
3. Inspect all connections.
4. Turn on vacuum source and test for any leakage.
5. Collect desired sample.

REMOVING SAMPLE

1. Remove tubing from Specimen Adaptor and cap port.
2. Remove Specimen Adaptor using caution.

DISPOSAL

- Dispose of according to your facility's standard operating procedures for biohazard disposal and both federal and local disposal regulations.

TROUBLESHOOTING

- Loss of Vacuum: Check that vacuum is on, canister is properly sealed, and that all connections are tight and tubing is not kinked. If loss continues, replace canister.

CAUTION

- Intended User: To be used by professionally trained healthcare personnel.
-  **Single Use Only. Do not attempt to clean, sterilize or reuse specimen adaptor.** Possible consequence of reuse include exposure to pathogens.
- **Do not** exceed vacuum level of 25" (640mm) Hg (84.7 kPa).
-  Contents are considered potentially hazardous. Use appropriate PPE and handle accordingly.
-  Store in a dark place. Long term exposure to light may compromise product performance and result in breakage during use.

Any serious incident related to the use of this product, should be reported to both the manufacturer and the health authority/competent authority where the product is established.

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Report any product malfunctions or complaints to Bemis Health Care at Customer.Service@BemisHealthCare.com

 Bemis Manufacturing Co.

**BEMIS**
HEALTH CARE™

 EMERGO EUROPE

CE

 Responsible Person: Bemis Limited