

Plastics Reference

	Suitable for Microwaves	Suitable for Autoclaving	Transparency	Max. Operating Temp (° C)	Brittleness Temp (° C)	Density (g/cm ³)
ETFE	Yes	Yes	Translucent	150	-100	1.70
FEP	Yes	Yes	Translucent	205	-100	2.15
HDPE	Yes	No	Translucent	105	-50	0.95
LDPE	Yes	No	Translucent	80	-50	0.92
PC	Yes	Yes*	Transparent	125	-130	1.20
PFA	Yes	Yes	Translucent	260	-200	2.17
PMMA	No	No	Transparent	65-95	-50	1.18
PMP	Yes	Yes	Transparent	150	0	0.83
POM	No	Yes*	Opaque	130	-40	1.42
PP	Yes	Yes	Translucent	125	0	0.90
PS	No	No	Transparent	70	-20	1.05
PTFE	Yes	Yes	Opaque	260	-200	2.17
SAN	No	No	Transparent	70	-40	1.03

* Frequent autoclaving may reduce mechanical stability. "Autoclaving" refers to steam sterilization at 121°C (250°F) at 15 PSIG for 20 minutes

+ Excellent chemical resistance	o Good to limited resistance	- Poor chemical resistance
Continuous exposure to the substance does not cause damage within 30 days. The plastic may remain resistant for years.	Continuous exposure to the substance causes minor damage, some of which is reversible, within 7-30 days (e.g., swelling, softening, decrease of mechanical strength, discoloration).	Not suitable for continuous exposure to the substance. Immediate damage may occur (loss of mechanical strength, deformation, discoloration, cracking, dissolution).

Chemical resistance of plastics to classes of substances at 20 °C

	SAN	PC	POM	PMP	LDPE	HDPE	PP	ETFE	PFA	FEP	PTFE	PMMA	PS
Alcohols, aliphatic	+	+	+	+	+	+	+	+	+	+	+	-	+
Ethers	-	-	+	-	o	o	o	+	+	+	+	-	-
Aldehydes	-	o	o	o	o	+	+	+	+	+	+	o	-
Esters	-	-	-	o	o	o	o	+	+	+	+	o	-
Hydrocarbons, aliphatic	-	o	+	o	o	+	+	+	+	+	+	+	-
Hydrocarbons, aromatic	-	-	+	-	o	+	o	+	+	+	+	-	-
Hydrocarbons, halogenated	-	-	+	-	o	o	o	+	+	+	+	-	-
Ketones	-	-	+	o	o	o	o	o	+	+	+	-	-
Alkalis	+	-	+	+	+	+	+	+	+	+	+	+	+
Acids, strong or concentrated	-	-	-	+	+	+	+	+	+	+	+	-	o
Acids, weak or diluted	o	o	-	+	+	+	+	+	+	+	+	-	o
Oxidizing acids, oxidizing agents	-	-	-	-	-	-	-	+	+	+	+	-	-

Abbreviations of the described plastics to DIN 7728

SAN	Styrene-acrylonitrile copolymer	PMP	Polymethylpentene
PMMA	Polymethyl methacrylate	ETFE	Ethylenetetrafluoro ethylene copolymer
PC	Polycarbonate	FEP	Perfluoroethylene-propylene copolymer
POM	Polyoxymethylene	PFA	Perfluoroalkoxy copolymer
LDPE	Low-density Polyethylene	PS	Polystyrene
HDPE	High-density Polyethylene	PTFE	Polytetrafluoroethylene
PP	Polypropylene		

The recommendations listed here are based on technical literature and information provided by the manufacturers of raw materials. They were prepared carefully and are intended to inform and advise. However, they cannot replace suitability testing performed by the user under actual working conditions.



Food-safe products according to EC Directive No. 10/2011