

# champion.

## OWNER'S OPERATING AND MAINTENANCE MANUAL



### Siesta Key

SLEEPER CHAIR



## **IMPORTANT**

**READ AND FAMILIARIZE YOURSELF WITH ALL INSTRUCTIONS BEFORE USING THIS PRODUCT!**

Do NOT install, maintain or operate this equipment without reading and following this manual, otherwise injury and/or damage may result.

IF YOU HAVE ANY QUESTIONS OR CONCERNS, PLEASE CONTACT CHAMPION MANUFACTURING, INC .

Champion assumes no responsibility for damage or injury caused by improper assembly, installation, use, or maintenance of these products.

No part of this manual may be duplicated in any form without the prior consent of Champion Manufacturing, Inc . Unauthorized duplication/distribution of these materials may result in civil prosecution to the maximum extent allowed by law.

The information contained in this manual is subject to change without notice.

**SAVE THESE INSTRUCTIONS FOR FUTURE REFERENCE!**

These instructions are available online at no charge.  
Visit [championchair.com/product-manual/](http://championchair.com/product-manual/) to download.

## **SYMBOLS**



**FOLLOW INSTRUCTIONS**



**GENERAL WARNING**



**PINCH-POINT WARNING**

# TABLE OF CONTENTS

Important Notes about this Manual.....	2
Transportation, Storage, Handling & Disposal.....	4
Preparation: Before You Begin.....	4
Cleaning Precautions .....	5
Service information .....	6
Determining the Problem .....	6
Serial Number .....	6
Parts Identification.....	6
Information Specific To The Siesta Key.....	
Operating Instructions.....	7
Upright To Chaise .....	7
Upright To Lay-Down.....	8
Schematics	
Upholstered Assemblies, Frame Base & Casters.....	9
Back Assembly .....	10
Seat Base Assembly .....	11
Specifications.....	12
Service & Contact Information .....	13

## TRANSPORTATION, STORAGE, HANDLING & DISPOSAL

- The product should be transported in factory packaging, inside an appropriate medium for the destination, i.e., air/sea cargo containers.
- The product should be stored in an environment that will inhibit rust/mold formation.
- Champion Manufacturing, Inc recommends not leaving the product in the factory packing in excess of three (3) months.
- The product should always be handled in a manner consistent with the user instructions and in a manner to prevent contamination after each use.
- The product has many recyclable components and to the extent practical, all effort should be used to recycle responsibly. Otherwise the product components shall be disposed of in accordance with local statutes.

## PREPARATION: BEFORE YOU BEGIN

1. Carefully examine your product for any damage. Be sure to inspect all components.  
**IF DAMAGE IS EVIDENT, CONTACT FREIGHT CARRIER OR CHAMPION IMMEDIATELY.**
2. Remove all packaging material and any hardware that was secured for shipping.
3. Carefully remove all components and any included tools and/or parts from the carton.
4. You may need to cut packaging materials with a box cutter or scissors to access the product. Use **CAUTION** to avoid personal injury or damage to the product.
5. Save all boxes and packaging material until **AFTER** you have assembled your product and have verified that all components are functioning properly. These material are required if it becomes necessary to return the product.
6. **DO NOT** install, maintain or operate this equipment without reading and following this manual otherwise injury and/or damage may result. **IF YOU HAVE QUESTIONS OR CONCERNS PLEASE CONTACT Champion Manufacturing, Inc .**

**CHAMPION ASSUMES NO RESPONSIBILITY FOR DAMAGE OR INJURY CAUSED BY IMPROPER ASSEMBLY, INSTALLATION, USE, OR MAINTENANCE OF THESE PRODUCTS.**

## GENERAL CLEANING PRECAUTIONS

### **WARNING** IMPORTANT: PLEASE READ

#### WARNINGS

- Treat all stains immediately to prevent it becoming a permanent stain.
- Several light applications of a cleaning agent are less harmful than a single concentrated application to remove a stain.
- **NEVER** mix ammonia, or a cleaner with ammonia, with bleach as dangerous compounds may result.
- When solvent type cleaners are being used, care should be exercised. **KEEP AWAY** from fire or flame and use in a well-ventilated area.
- Some fabric dyes, including indigo dyes commonly used in denim jeans, may migrate to the surface of the material and have the ability to create a permanent stain. This is increased by humidity and temperature and may be irreversible. To reduce the likelihood of a problem from this or any stain, clean as soon as noticed.
- Indiscriminate use of strong solvents can damage or discolor vinyl.
- Do not use abrasive household cleansers or steel wool.
- Residue left by cleaning agents will degrade and shorten the life and/or affect the product's appearance if not removed.
- Do not use steam cleaning or hot water to clean.
- Do not use an iodine-based solution or the upholstery will become stained.

## SERVICE INFORMATION

The mission of the customer service department is to get your chair up and running as quickly as possible. It is critical that the customer service department know what product you have, and exactly what is wrong with the product. If you have questions or problems, you should never hesitate to call for assistance: 800-237-3377.

The most timely and cost-effective way for your chair to be repaired is for the Champion's customer service department to work with your maintenance department or equipment technician.

## DETERMINING THE PROBLEM

What is wrong with the chair should be determined by troubleshooting. The customer service department will assist you with this by asking you questions about the chair function.

## SERIAL NUMBER

The chair serial number identifies the precise configuration of your chair; this is critical to receiving correct components and instructions. This number is required to process your request.

**On the Siesta Key sleeper:** to find your serial number, go to back under lower left flap, separate the hook and loop fastener to see the serial number.

## PARTS IDENTIFICATION

To identify worn or damaged components please refer to appropriate product schematics.

To obtain repair part numbers refer to the parts listing key using the schematics page and item number and please call Champion Manufacturing at 800-237-3377 for assistance.

Parts orders may be placed by using the convenient fax order form in this manual or by calling Champion's customer service at 800-237-3377.

When placing an order by phone you will be asked to provide:

- your name and facility name
- your phone, fax number, email address
- the serial number of your product, and
- the nature of your problem

Having the above information available at the time you call will expedite the process. In order to provide prompt, accurate service it may be necessary to request further information about the chair function to accurately define the problem.

# OPERATING INSTRUCTIONS

## CREATING UPRIGHT TO CHAISE POSITION






### WARNING

This chair has moving parts that may create pinch points. Keep your hands clear of the mechanism when you move the chair from one position to another.

1. Lock rear casters.
2. Raise seat cushions. Use the lift strap to pull the seat base up and out. Place seat cushion behind the seat base.
3. To return to upright, reverse the process. Lift the seat cushion, then use the strap to lift the seat base up and in, to the returned position. Place the seat cushion on top of the seat base.
4. To go to the lay down position, after extending the seat base, see the illustration on the next page.

*This illustration shows how to place a Siesta Key in chaise position. Also use the illustration to identify basic chair parts.*



### CHAIR PARTS

A	Back Pillow
B	Seat Cushion
C	Mid-Section
D	Lift Strap
E	Front Non-Locking Caster
F	Rear Locking Caster
G	Upholstered Arm

# OPERATING INSTRUCTIONS – CONTINUED



## CREATING UPRIGHT TO LAY-DOWN POSITION



### WARNING

This chair has moving parts that may create pinch points. Keep your hands clear of the mechanism when you move the chair from one position to another.


1. Lock rear casters.
2. Raise seat cushions. Use the lift strap to pull the seat base up and out. Move to the back of the chair. Pull outward on the back kick-out to release the chair back. Allow the chair back to move back until it is supported by the kick-out.
3. Finish by placing the seat cushion between the extended seat base and back.
4. To return to upright, lift the seat cushion, then use the strap to lift the seat base up and back in, to the returned position. Place the seat cushion on top of the seat base. Finish by moving to the back of the chair and lifting until the kick-out release snaps back into place.



*This illustration shows how to place a Siesta Key in the lay-down position. Also use the illustration to identify basic chair parts.*

### CHAIR PARTS

A	Back Pillow
B	Seat Cushion
C	Mid-Section
D	Lift Strap
E	Front Non-Locking Caster
F	Rear Locking Caster
G	Upholstered Arm
H	Back Kick-out



# SCHEMATICS


## UPHOLSTERED ASSEMBLIES, FRAME BASE AND CASTERS



**WARNING**

This chair has moving parts that may create pinch points. Keep your hands clear of the mechanism when you move the chair from one position to another.

Use the illustration to identify basic chair parts.



### CHAIR PARTS

A	Back Pillow Assembly
B	Arm Asm.; LH, Shown
C	Rear Locking Caster
D	Front Swivel Caster
E	Frame Base
F	Base Insert
G	Seat Base Pull-Out Asm.
H	Seat Cushion Assembly
H1	Seat Cushion Cover Only
H2	Seat Foam Only (Not shown)

## SCHEMATICS – CONTINUED


### BACK ASSEMBLY



#### WARNING

This chair has moving parts that may create pinch points. Keep your hands clear of the mechanism when you move the chair from one position to another.

Use the illustration to identify basic chair parts.



#### CHAIR PARTS

A	Back Pillow Assembly
A1	Back Pillow Cover Only (Not shown)
A2	Back Pillow Foam Only (Not shown)
B	Back Frame Assembly
C	Foam: Top Strip
D	Outside Back Upholstered Asm.
E	Back Kick-Out Uph. Cover
F	Back Fabric Panel

# SCHEMATICS – CONTINUED


## SEAT BASE ASSEMBLY



**WARNING**

This chair has moving parts that may create pinch points. Keep your hands clear of the mechanism when you move the chair from one position to another.

Use the illustration to identify basic chair parts.




### CHAIR PARTS

A	Seat Base Assembly
A1	Seat Base Cover (Not shown)
A2	Seat Base Foam (Not shown)
B	Seat Base Frame
C	Pull-Out Lever
D	Pull Strap
E	Seat Base Face Board Asm.
F	Glide, Front Leg
G	Leg Set, Frt. Face Board

# IMPORTANT: ALWAYS REFER TO THE WEIGHT CAPACITY LABEL ON YOUR CHAIR.

In all cases the labeling on the chair at the time of delivery indicates the correct rating for your chair – weight rating should not be exceeded!

## SPECIFICATIONS



### Siesta Key

Weight capacity	300 lb. (136,078 kg)
Overall height	42" (106.68 cm)
Overall width	32.5" (82.55 cm)
Sleep length	74" (187.96 cm)
Seat height/width	18.75"/25" (45.72 cm – 63.5 cm)
Seat depth	21" (53.34 cm)
Arm width	2.5" (6.35 cm)
Optional wood armrest width	3" (7.62 cm)
Total number of positions	3

**IMPORTANT: Maximum Weight Capacity should be evenly distributed.**