

Instructions for Users

Model: PULSAR PRO LOW

PULSAR PRO LOW

PULSAR PRO LOW is intended for the support of mild to moderate ankle distortions, Tibialis Tendinitis, Osteoarthritis of ankle and foot, Chronic ankle instability, Tendinopathy of foot and ankle; Prophylactic and therapeutic use in sports activities. Rupture - tendon of the foot and ankle, non-traumatic. Acquired flat-foot or pronation of the foot/ankle.



1. Open all the straps and slide the injured ankle into the brace.



2. Put the foot of the injured ankle upon the bottom foot plate and secure the strap on the foot.



3. Locate the top strap and feed through the D-Ring on the opposite side and attach to desired support. The length of this strap is adjustable by shortening or lengthening the strap on the underside. DO NOT OVERTIGHTEN.



CARE AND CLEANING INSTRUCTIONS

Hand wash with a mild soap and let air dry.

PRO ANKLE BRACE SIZE CHART

Size	Women's Shoe Size	Men's Shoe Size
S/M Left	Up to Size 10	Up to Size 9
S/M Right	Up to Size 10	Up to Size 9
L/XL Left	10.5+	9.5+
L/XL Right	10.5+	9.5+