

Forearm Crutches

Muletas de Lujo de Antebrazo

Items 6350A, 6350T, 6350Y



SAFETY GUIDELINES - PLEASE READ BEFORE USE

1. This product should not be used without proper instruction from a healthcare professional.
2. This product has a maximum user weight limitation of 300 lb. - EVENLY DISTRIBUTED.
3. GF Health Products, Inc. assumes no responsibility for any damage or injury caused by improper installation, assembly or use of this product.
4. If components are damaged or missing, contact your distributor immediately. DO NOT use substitute parts.
5. Before use, ensure tip is properly seated and straight. Check rubber tip for rips, tears, cracks or wear. If any of these conditions exist, replace rubber tip immediately.
6. Always use caution when using this product on wet or slippery surfaces. DO NOT use crutches on stairs or steep inclines.
7. Inspect product weekly for proper function and wear. If any component is not in proper working order or appears worn, contact your Lumex dealer for repair and/or appropriate replacement parts.
8. Before use, ensure that push buttons protrude fully through adjustment holes.
9. Make sure locking collet is securely tightened at all times.

GENERAL WARNING

- **DO NOT** use this product without first completely reading and understanding these instructions and any additional instructional material such as owner's manuals, service manuals or instruction sheets supplied with this product. If you are unable to understand the warnings, cautions or instructions, contact a healthcare professional, dealer or technical personnel before attempting to use this equipment - otherwise, injury or damage may occur.
- If components are damaged or missing, contact your GF Health Products, Inc. distributor immediately. **DO NOT** substitute parts.

GENERAL WARNING (CONT.)

- Before use, ensure that all assembled parts are secure and all hardware is tightened.
- Care should be taken to ensure that all height adjustments are secure.
- **ALWAYS** observe the user weight limitation on the hangtag of this product.

 **WARNING:** Cancer and Reproductive Harm-
www.p65warnings.ca.gov.

INSTALLATION WARNINGS

- GF Health Products, Inc. assumes no responsibility for any damage or injury caused by improper installation, assembly or use of this product.
- Inspect product weekly or on a set routine preventative maintenance schedule for normal wear and tear. If any component is not in proper working order or appears worn, contact your Lumex dealer for repair and/or appropriate replacement parts.
- **ALWAYS** test to see that the product is properly secured in open position **BEFORE** use and ensure all hardware is securely tightened.

MAINTENANCE

- Wipe product weekly with a damp cloth.
- Grips may be cleaned with mild soap or detergent. **DO NOT** use organic solvents, such as acetone, which may damage grip.
- Never use organic solvents, abrasive cleaners, or scouring pads on any part of the product.
- Check hardware and positioning weekly for tightness.
- Replace crutch tips as needed.
- To guarantee the best service from your product, always specify Lumex replacement parts, available from your GF Health Products, Inc. distributor.

INSTRUCTIONS FOR USE

Please Read all Instructions Before Use.

The chart below can be used as a guide to select height-appropriate crutches. Please consult a healthcare professional for assistance with proper height adjustment.

Item Number	Size	User Height	Adjustment Range / Handgrip to Floor
6350A	Adult	5'0" - 6'0" (1.5 - 1.8m)	28" - 37" (71 - 94 cm)
6350T	Tall	5'5" - 6'4" (1.65 - 1.93m)	33" - 42" (83.8 - 106.7 cm)
6350Y	Youth	4'10" - 5'8" (1.47 - 1.72m)	22" - 31" (55.9 - 78.7 cm)

ADJUST CRUTCH HEIGHT

1. Loosen lower locking collet by twisting counter-clockwise.
2. Adjust crutch height: Depress lower push button and slide lower section upward or downward to desired height.

NOTE: *An audible "click" will be heard when push button securely protrudes through the adjustment hole.*

3. Secure lower locking collet by turning clockwise until tight.

ADJUST CUFF HEIGHT

1. Loosen upper locking collet by twisting counter-clockwise.
2. Adjust cuff height: Depress upper push button and slide upper section upward or downward to desired height.

NOTE: *An audible "click" will be heard when push button securely protrudes through the adjustment hole.*

3. Secure upper locking collet by turning clockwise until tight.
4. Arm cuffs are adjustable; adjust both cuffs for proper fit.

LIMITED WARRANTY

GF Health Products, Inc. ("GF") warrants this product to be free from defects in materials and workmanship to the original retail purchaser or original lessee only ("Customer") for:

Frame: three (3) years

Tips, Wheels/Casters, and Handgrips: three (3) months

The warranty does not cover: 1) other non-durable components; 2) defects, damage or other conditions caused by misuse, abuse, negligence, alteration, accident, freight damage, tampering or failure to seek and obtain repair or replacement in a timely manner; and 3) products not installed, used, or properly cleaned and maintained in accordance with the applicable GF instructions. Product provided to a lessee must be inspected prior to delivery and must be maintained by the Dealer. The lessee must be properly trained on the operation and safety of this product prior to use.

GF will, at its option, repair or replace the defective product. Warranty issues must: 1) be reported to the Dealer from whom the Customer purchased the product within the applicable warranty period and 2) upon examination by GF or its authorized representative, prove to be a warranty item. GF is not responsible for freight or labor charges associated with the warranted item. If there is no Dealer, you must contact GF directly by calling 770-368-4700, faxing 770-368-4932 or e-mailing cs@grahamfield.com. Specific directions will be provided.

TO THE GREATEST EXTENT PERMITTED BY LAW, THIS WARRANTY IS GF'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. EXCEPT IN CASES WHERE IMPLIED WARRANTIES MAY NOT BE VALIDLY WAIVED, GF MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. GF'S TOTAL LIABILITY FOR THIS PRODUCT IS LIMITED TO THE COST OF THE PRODUCT(S) GIVING RISE TO THE CLAIM. IN NO EVENT WILL GF BE LIABLE FOR ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES. CERTAIN LAWS MAY CONFER ADDITIONAL RIGHTS REGARDING WARRANTIES FOR SPECIFIC PRODUCTS AND IN THOSE INSTANCES GF'S LIABILITY AND THE LIABILITY OF GF'S SUPPLIERS, SHALL BE LIMITED TO THE FULLEST EXTENT PERMITTED BY LAW. This warranty gives you specific legal rights. You may also have other rights, which vary from state to state or country to country.

Cuff Assembly
Ensamblaje del puño

Upper Locking Collet
Cerradura Superior

Upper Push Button
Botón de Apretar Superior

Handgrip
Manga

Lower Push Button
Botón de Apretar Inferior

Lower Locking Collet
Cerradura Inferior

Tip
Punta



Features

- Soft, cushion grip reduces hand stress and fatigue
- Leg and forearm sections can be adjusted independently for a custom fit
- Vinyl-coated arm cuffs can be shaped to the user's arm for extra comfort
- Extra-large, non-skid crutch tips provide added traction and stability
- Easy push-button height adjustment with locking collet for added security
- Grip-to-floor height adjusts in 1" increments
- Arm cuff to grip height adjusts in 1" increments
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED

