



PL400 / PL600 Digital Patient Lift Scale

Patient Lift Scale is a completely self-contained instrument that determines the weight of a person suspended from a portable patient lift mechanism. It operates on only one (1) 9-volt Alkaline battery (not included). Battery power is saved due to the automatic shut-off feature. The scale also has a membrane keypad with only three (3) keys: ON/ZERO, OFF and LB/KG. It's so easy to setup and operate; the procedure is quickly learned.



PL400/600 Specifications

Capacity:	PL400 400 lb x .2 lb / 180 kg x .1 kg PL600 600 lb x .2 lb / 270 kg x .1 kg
Power Source:	One (1) nine-volt battery (not included)
Keypad Controls:	on/zero, off, and lb/kg display select
Display:	.75" (19mm) high-contrast LCD with annunciator for lb/kg, zero, low battery and display lock
Size:	4.125" W x 2.125" D x 4" H / 104mm W x 54mm D x 102mm H (center of eye to center of eye: 5" overall - 127 mm)