

Distributed by: Universal Medical Inc.

www.universalmedicalinc.com



FHD-133-II / FHD-144-II

Geriatric/Bariatric Wheelchair Scale



The battery operated models FHD133II and FHD144II wheelchair scales have a nonskid, rubber platform and ramp. The large platform accommodates walk-on patients, wheelchairs, or straight-legged chairs. This large capacity scale can be mounted flush with the floor if desired. The remote digital weight display includes a .7" (18mm) high LCD readout with keypad tare and is suitable for desk or wall mount.

FHD133II / FHD144II Specifications

Capacity:	1,000 lb x .2 lb/450 kg x .1 kg
Platform Sizes:	FHD133II 36" x 36" (914mm x 914mm) FHD144II 48" x 48" (1.21m x 1.21m)

Platform Material:	Smooth steel covered with black rubberized mat
Shipping Weights:	FHD133II 315 lb (143 kg) FHD144II 460 lb (209 kg)
Ramp:	28.5" / 72 cm ramp (one included)
Display:	.7" (18mm) high contrast LCD
Cable Length:	15' / 4.5 m
Color:	Blue with black mat and ramp