



**[HM2002-38EXT Solo Spot  
C Clamp Exam Light](#)**

**[HM2005-37EXT Solo Spot  
Wall Mounted Exam Light](#)**

**[HM2006-38EXT Solo Spot  
Mobile Exam Light](#)**

**[HM2017-44 Solo Spot  
Ceiling Mount Exam Light](#)**

The SoloSpot is ideal for applications where a focused, spot light is required. The specially designed lamp head offers cool light and a focused, intense spot. This light is perfect for direct attachment to medical equipment such as patient tables and chairs. Ideal for OB-GYN, dental, and surgical applications.

**SOLOSPOT™ FEATURES**

>>> Unique, double filtered lens and color correction filter with lampshade extension help keep light cool on skin and illuminate with natural skin tones to improve tissue clarity

- >>> Excellent maneuverability allows for easy adjustment of light
- >>> Intense 20W Halogen light source (average 2,000 hours)
- >>> 3.5-inch diameter spot pattern @ 20-inches
- >>> Available in wall, ceiling, c-clamp and mobile configurations



Mobile

**TECHNICAL DATA**

**Power supply:** 120V 60 Hz AC or 230-240V 50 Hz AC

**Light intensity:** Up to 3,000 Footcandles/32,400 Lux

**Light source:** Halogen bulb – 12V/20W

**Light color temp:** 4,000 Kelvin

**Color:** White – RAL 9002

| Configurations: | Wall  | Ceiling | C-clamp | Mobile |
|-----------------|-------|---------|---------|--------|
| Weight:         | 6 lbs | 16 lbs  | 6 lbs   | 15 lbs |

|         |                               |
|---------|-------------------------------|
| Wall    | Part No.: <b>HM2005-37EXT</b> |
| C-clamp | Part No.: <b>HM2002-38EXT</b> |
| Mobile  | Part No.: <b>HM2006-38EXT</b> |
| Ceiling | Part No.: <b>HM2017-44EXT</b> |



C-Clamp



Wall