



7156-70: 2 Pegs. 2 mounting holes denoted with dots.  
 1. located 1" from the top centered at 3" from either side and 2. located 1" from the bottom centered at 3" from either side.




7156-80: 4 Pegs. 2 mounting holes denoted with dots. : 1. located 1" from the top, 3" from the left. 2. located 1" from the top, 3" from the right.



7156: 8 Pegs. 4 mounting holes (dots): 1 & 4. located 1" from the top, 3 1/8" from the left and right. 2 & 3. located 1" from the top, 9 1/8" from the left and right



7156-35: 14 Pegs. 5 mounting holes (dots): 1 & 5 located 1 1/8" from the top, 2 1/16" from the left and right. 2 & 4 located 1 1/8" from the top, 6 1/16" from the left and right. 3 located 1 1/8" from the top, 18" from either side.



7156-46: 18 Pegs. 5 mounting holes (dots): 1 & 5 located 1 1/8" from the top, 1 1/8" from the left and right. 2 & 4 located 1 1/8" from the top, 17 1/8" from the left and right. 3 located 1 1/8" from the top, 25" from either side.

