

SAFETY DATA SHEET

according to US Regulation 29 CFR 1910.1200 and the Canadian HPA

SONO Hydrogen Peroxide Disinfectant

Version 1.4

Revision Date 2020.04.20

Print Date 2020.07.24

SECTION 1. IDENTIFICATION

Product name : SONO Hydrogen Peroxide Disinfectant

PMRA Registration number : 6836-385-89018

Manufacturer or supplier's details

Company : Advanced Ultrasound Solutions, Inc.
1012 W. Beverly Blvd, #334
Montebello, CA 90640
USA

Telephone : 1-855-879-4737

E-mail address : info@sonowipes.com

Emergency telephone number : For incidents only (spill, leak, fire, exposure, or accident), call CHEMTREC at
1-800-424-9300 (inside North America) [CCN 864796]
1-703-741-5970 (outside North America) [CCN 864796]

+41 61 313 94 94 (24h)

SECTION 2. HAZARDS IDENTIFICATION

GHS Classification

Not a hazardous substance or mixture according to US Regulation 29 CFR 1910.1200 and the Canadian HPA.

GHS label elements

Based on available data, the classification criteria are not met.

Handle in accordance with good industrial hygiene and safety practice. Not a hazardous substance or mixture according to US Regulation 29 CFR 1910.1200 and the Canadian HPA.

Precautionary statements : **Prevention:**
P261 Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P262 Do not get in eyes, on skin, or on clothing.
P270 Do not eat, drink or smoke when using this product.
P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.
Response:
P314 Get medical advice/ attention if you feel unwell.
Storage:
P402 + P404 Store in a dry place. Store in a closed container.
P410 + P403 Protect from sunlight. Store in a well-ventilated place.
Disposal:
P501 Dispose of contents/container in accordance with local regulation.

Other hazards

None known.

SONO Hydrogen Peroxide Disinfectant

SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS

Hazardous components

Chemical name / Synonyms	CAS-No.	Concentration (% w/w)
Methanesulphonic acid	75-75-2	1 - 2.5
Hydrogen peroxide	7722-84-1	1 - 2.5
Orthophosphoric acid	7664-38-2	0.5 - 1

SECTION 4. FIRST AID MEASURES

If inhaled	: Move to fresh air. If unconscious, place in recovery position and seek medical advice.
In case of skin contact	: Wash with water and soap as a precaution. If on clothes, remove clothes. In the case of skin irritation or allergic reactions see a physician.
In case of eye contact	: Flush eyes with water as a precaution. Remove contact lenses. Keep eye wide open while rinsing. Protect unharmed eye. If eye irritation persists, consult a specialist.
If swallowed	: Clean mouth with water and drink afterwards plenty of water. Never give anything by mouth to an unconscious person. If symptoms persist, call a physician.
Most important symptoms and effects, both acute and delayed	: No information available.
Notes to physician	: Treat symptomatically.

SECTION 5. FIREFIGHTING MEASURES

Suitable extinguishing media	: Water spray Alcohol-resistant foam Dry chemical
Unsuitable extinguishing media	: High volume water jet
Specific hazards during firefighting	: Heating or fire can release toxic gas.
Further information	: Use water spray to cool unopened containers.
Special protective equipment for firefighters	: In the event of fire, wear self-contained breathing apparatus. Use personal protective equipment.

SONO Hydrogen Peroxide Disinfectant

SECTION 6. ACCIDENTAL RELEASE MEASURES

- Personal precautions, protective equipment and emergency procedures : Use personal protective equipment.
- Environmental precautions : Try to prevent the material from entering drains or water courses.
- Methods and materials for containment and cleaning up : Neutralize with chalk, alkali solution or ammonia. Contain spillage, and then collect with non-combustible absorbent material, (e.g. sand, earth, diatomaceous earth, vermiculite) and place in container for disposal according to local / national regulations (see section 13).

SECTION 7. HANDLING AND STORAGE

- Advice on protection against fire and explosion : Normal measures for preventive fire protection.
- Advice on safe handling : No special precautions required.
- Conditions for safe storage : Keep container tightly closed. Keep in a well-ventilated place. Electrical installations / working materials must comply with the technological safety standards. To maintain product quality, do not store in heat or direct sunlight.
- Materials to avoid : No special restrictions on storage with other products.
- Further information on storage stability : No decomposition if stored and applied as directed.

SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Components with workplace control parameters

Components	CAS-No.	Value type (Form of exposure)	Control parameters / Permissible concentration	Basis
Hydrogen peroxide	7722-84-1	TWA	1 ppm	ACGIH
		REL	1 ppm 1.4 mg/m ³	NIOSH/GUIDE
Orthophosphoric acid	7664-38-2	TWA	1 mg/m ³	ACGIH
		STEL	3 mg/m ³	ACGIH
		REL	1 mg/m ³	NIOSH/GUIDE
		STEL	3 mg/m ³	NIOSH/GUIDE

Personal protective equipment

- Respiratory protection : No personal respiratory protective equipment normally required.

SONO Hydrogen Peroxide Disinfectant

Hand protection

Remarks	: Handle in accordance with good industrial hygiene and safety practice.
Eye protection	: Handle in accordance with good industrial hygiene and safety practice.
Skin and body protection	: Handle in accordance with good industrial hygiene and safety practice.
Hygiene measures	: Handle in accordance with good industrial hygiene and safety practice.

SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance	: liquid (77 °F / 25 °C)
Colour	: clear
Odour	: Detergent
Odour Threshold	: no data available
pH	: 1.07 (77 °F / 25 °C) (undiluted)
Melting point/freezing point	: no data available
Boiling point/boiling range	: no data available
Flash point	: > 212 °F / > 100 °C
Evaporation rate	: no data available
Flammability (solid, gas)	: no data available
Flammability (liquids)	: no data available
Self-ignition	: no data available
Upper explosion limit	: no data available
Lower explosion limit	: no data available
Vapour pressure	: no data available
Relative vapour density	: no data available
Relative density	: no data available
Density	: 0.0137 g/cm ³ (68 °F / 20 °C)
Bulk density	: no data available

SONO Hydrogen Peroxide Disinfectant

Water solubility	:	soluble (77 °F / 25 °C)
Partition coefficient: n-octanol/water	:	no data available
Auto-ignition temperature	:	no data available
Decomposition temperature	:	no data available
Viscosity, dynamic	:	no data available
Viscosity, kinematic	:	0.98 mm ² /s (77.2 °F / 25.1 °C)
Explosive properties	:	no data available
Oxidizing properties	:	no data available
Minimum ignition energy	:	no data available

SECTION 10. STABILITY AND REACTIVITY

Reactivity	:	No decomposition if stored and applied as directed.
Chemical stability	:	Stable under recommended storage conditions.
Possibility of hazardous reactions	:	Stable under recommended storage conditions. No hazards to be specially mentioned.
Conditions to avoid	:	Heat
Incompatible materials	:	Strong acids and strong bases Oxidizing agents
Hazardous decomposition products	:	No decomposition if used as directed.

SECTION 11. TOXICOLOGICAL INFORMATION

Acute toxicity

Acute oral toxicity	:	LD50 (Rat, female): > 5,000 mg/kg Method: US-EPA GLP: yes
Acute inhalation toxicity	:	LC50 (Rat, male and female): > 2.52 mg/l Exposure time: 4 h Test atmosphere: dust/mist Method: US-EPA GLP: yes Assessment: The substance or mixture has no acute inhalation toxicity
Acute dermal toxicity	:	LD50 (Rat, male and female): > 5,000 mg/kg Method: US-EPA GLP: yes

Skin corrosion/irritation

SONO Hydrogen Peroxide Disinfectant

Species: Rabbit
 Method: US-EPA
 Result: No skin irritation
 GLP: yes

Serious eye damage/eye irritation

Species: Rabbit
 Result: No eye irritation
 Method: US-EPA
 GLP: yes

Respiratory or skin sensitisation

Species: Guinea pig
 Method: US-EPA
 Result: Not a skin sensitizer.
 GLP: yes

Germ cell mutagenicity

Genotoxicity in vitro : Remarks: no data available

Carcinogenicity

Remarks: no data available

IARC

No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

OSHA

No component of this product present at levels greater than or equal to 0.1% is on OSHA's list of regulated carcinogens.

NTP

No component of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

ACGIH

Confirmed animal carcinogen with unknown relevance to humans
 Hydrogen peroxide 7722-84-1

Reproductive toxicity

Effects on fertility : Remarks: no data available

STOT - single exposure

Remarks: no data available

STOT - repeated exposure

Remarks: no data available

Aspiration toxicity

No aspiration toxicity classification

Further information

Remarks: no data available

SONO Hydrogen Peroxide Disinfectant

SECTION 12. ECOLOGICAL INFORMATION

Ecotoxicity

Toxicity to fish : Remarks: no data available

Persistence and degradability

no data available

Bioaccumulative potential

Bioaccumulation : Remarks: no data available

Components:

Methanesulphonic acid:

Partition coefficient: n-octanol/water : log Pow: -4.98

Mobility in soil

Distribution among environmental compartments : Remarks: no data available

Other adverse effects

Ozone-Depletion Potential : Regulation: US. EPA Clean Air Act (CAA) Section 602 Ozone-Depleting Substances (40 CFR 82, Subpt. A, App A & B)
Remarks: This product neither contains, nor was manufactured with a Class I or Class II ODS as defined by the U.S. Clean Air Act Section 602 (40 CFR 82, Subpt. A, App.A + B).

Additional ecological information : no data available

SECTION 13. DISPOSAL CONSIDERATIONS

Disposal methods

Waste from residues : Dispose of contents/container in accordance with local regulation.
Contact waste disposal services.

Contaminated packaging : Empty containers should be taken to an approved waste handling site for recycling or disposal.

SECTION 14. TRANSPORT INFORMATION

DOT : Not dangerous goods

UN number : Not applicable
Proper shipping name : Not applicable
Transport hazard class : Not applicable
Packing group : Not applicable

SONO Hydrogen Peroxide Disinfectant

TDG	:	Not dangerous goods
UN number	:	Not applicable
Proper shipping name	:	Not applicable
Transport hazard class	:	Not applicable
Packing group	:	Not applicable
IATA	:	Not dangerous goods
UN number	:	Not applicable
Proper shipping name	:	Not applicable
Transport hazard class	:	Not applicable
Packing group	:	Not applicable
IMDG	:	Not dangerous goods
UN number	:	Not applicable
Proper shipping name	:	Not applicable
Transport hazard class	:	Not applicable
Packing group	:	Not applicable
ADR	:	Not dangerous goods
UN number	:	Not applicable
Proper shipping name	:	Not applicable
Transport hazard class	:	Not applicable
Packing group	:	Not applicable
RID	:	Not dangerous goods
UN number	:	Not applicable
Proper shipping name	:	Not applicable
Transport hazard class	:	Not applicable
Packing group	:	Not applicable
Special precautions for user	:	none
Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code	:	Not applicable

SECTION 15. REGULATORY INFORMATION

This chemical is a pesticide product registered by the United States Environmental Protection Agency and is subject to certain labeling requirements under federal pesticide law. These requirements differ from the classification criteria and hazard information required for safety data sheets (SDS), and for workplace labels of non-pesticide chemicals.

EPA Registration number	:	6836-385-89018
Signal word	:	
Hazard statements	:	

This chemical is a pest control product registered by Health Canada Pest Management Regulatory Agency and is subject to certain labelling requirements under the Pest Control Products Act. These requirements

SONO Hydrogen Peroxide Disinfectant

differ from the classification criteria and hazard information required for GHS-consistent safety data sheets. Following is the hazard information required on the pest control product label:

PMRA Registration number : 6836-385-89018
 Hazard pictograms :
 Signal word :
 Hazard statements :

EPCRA - Emergency Planning and Community Right-to-Know Act

CERCLA Reportable Quantity

Components	CAS-No.	Component RQ (lbs)	Calculated product RQ (lbs)
Orthophosphoric acid	7664-38-2	5000	*

*: Calculated RQ exceeds reasonably attainable upper limit.

SARA 304 Extremely Hazardous Substances Reportable Quantity

Components	CAS-No.	Component RQ (lbs)	Calculated product RQ (lbs)
Hydrogen peroxide	7722-84-1	1000	*

*: Calculated RQ exceeds reasonably attainable upper limit.

SARA 311/312 Hazards

See above: SECTION 2. Hazard Identification-GHS Classification

SARA 302 Extremely Hazardous Substances Threshold Planning Quantity

Components	CAS-No.	Concentration
Hydrogen peroxide	7722-84-1	>= 1 - < 5 %

SARA 313

This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

Clean Air Act

This product neither contains, nor was manufactured with a Class I or Class II ODS as defined by the U.S. Clean Air Act Section 602 (40 CFR 82, Subpt. A, App.A + B).

This product does not contain any hazardous air pollutants (HAP), as defined by the U.S. Clean Air Act Section 112 (40 CFR 61).

This product does not contain any chemicals listed under the U.S. Clean Air Act Section 112(r) for Accidental Release Prevention (40 CFR 68.130, Subpart F).

This product does not contain any chemicals listed under the U.S. Clean Air Act Section 111 SOCM I Intermediate or Final VOC's (40 CFR 60.489).

This product does not contain any VOC exemptions listed under the U.S. Clean Air Act Section 450.

SONO Hydrogen Peroxide Disinfectant

Clean Water Act

The following Hazardous Chemicals are listed under the U.S. CleanWater Act, Section 311, Table 117.3:

Components	CAS-No.	Component RQ (lbs)
Orthophosphoric acid	7664-38-2	5000

The following Hazardous Substances are listed under the U.S. CleanWater Act, Section 311, Table 116.4A:

Components	CAS-No.	Concentration
Orthophosphoric acid	7664-38-2	>= 0.1 - < 1 %

This product does not contain any toxic pollutants listed under the U.S. Clean Water Act Section 307

US State Regulations

Massachusetts Right To Know

Components	CAS-No.
Hydrogen peroxide	7722-84-1

Pennsylvania Right To Know

Components	CAS-No.
Water	7732-18-5

California Prop. 65

This product does not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

Canadian lists

NPRI

Canadian National Pollutant Release Inventory (NPRI): No component is listed on NPRI.

SECTION 16. OTHER INFORMATION

Full text of other abbreviations

ACGIH : US. ACGIH Threshold Limit Values
 NIOSH/GUIDE : US. NIOSH: Pocket Guide to Chemical Hazards, as amended

AICS - Australian Inventory of Chemical Substances; ASTM - American Society for the Testing of Materials; bw - Body weight; CERCLA - Comprehensive Environmental Response, Compensation, and Liability Act; CMR - Carcinogen, Mutagen or Reproductive Toxicant; DIN - Standard of the German Institute for Standardisation; DOT - Department of Transportation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; EHS - Extremely Hazardous Substance; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; HMIS - Hazardous Materials Identification System; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation

SONO Hydrogen Peroxide Disinfectant

Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; MSHA - Mine Safety and Health Administration; n.o.s. - Not Otherwise Specified; NFPA - National Fire Protection Association; NO(A)EC - No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NTP - National Toxicology Program; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; RCRA - Resource Conservation and Recovery Act; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; RQ - Reportable Quantity; SADT - Self-Accelerating Decomposition Temperature; SARA - Superfund Amendments and Reauthorization Act; SDS - Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative

Revision Date : 2020.04.20

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

Date format : yyyy/mm/dd

US / EN