

OEM Equipment Compatibility Chart by Imaging Modality

Manufacturer	Model	Attachment Method	Modality	Table Dimensions
GE	Brightspeed	4mm cording	CT	
GE	Lightspeed	4mm cording	CT	
GE	LX	Scan-Track™	CT	
GE	VCT	4mm or 5mm cording	CT	
GE	Prospeed Plus	Scan-Track™	CT	13" x 79"
GE	Hilight Advantage/9800	Scan-Track™	CT	15" x 78"
GE	Hilight Advantage/9800	Scan-Track™	CT	16" x 73"
GE	VCT1700	Scan-Track™	CT	17" x 83"
GE	VCT2000	Scan-Track™	CT	17" x 102"
GE	Innova 2100	5mm cording	Cath	
GE	Optima	Scan-Track™	CT	17" x 102"
GE	Optima 660	Scan-Track™	CT	17" x 79"
GE	Lite	Scan-Track™	CT	18" x 86"
GE	Discovery	Scan-Track™	CT	18" x 105"
Philips	Tomoscan	Scan-Track™	CT	10" x 69"
Philips	Brightview	Scan-Track™	CT	
Philips	Brilliance	5mm cording	CT	
Philips	Ingenuity	Scan-Track™	CT	
Siemens	128S Plus	4mm cording	CT	
Siemens	Sensation	5mm cording	CT	
Siemens	Somatom Sensation 40	5mm cording	CT	94.5" x 17.75"
Siemens	Somatom Definition AS	5mm cording	CT	18" x 94"
Siemens	Somatom Drive	5mm cording	CT	16" x 82"
Siemens	Somatom Force	5mm cording	CT	
Siemens	Somatom AR Star	5mm cording	CT	17" x 77"
Siemens	Somatom Emotion	Scan-Track™	CT	18" x 86"
Siemens	Somatom S Plus	5mm cording	CT	17" x 87"
Canon/Toshiba	Aquilion One	Scan-Track™	CT	96.5" x 19"
Canon/Toshiba	Asteion	Scan-Track™	CT	74" x 18"
GE	Hawkeye Millennium	Scan-Track™	NM	14" x 79"
GE	Infinia	Scan-Track™	NM	14" x 79"
GE	Ventri	Scan-Track™	NM	16" x 64"
GE	Myocite	Scan-Track™	NM	
Siemens	Symbia S	Scan-Track™	NM	18" x 80"
Siemens	Symbia T	Scan-Track™	NM	18" x 80"
Siemens	Symbia Evo	Scan-Track™	NM	