

Real Time Instrument Brand	Model Number/Name	Compatible Accuris qMax™ Reagents	
Applied BioSystems®	7000, 7300, 7700	qMax Green High Rox	qMax Probe High Rox
	7900, 7900HT, 7900HT Fast		
	StepOne™, StepOne™ Plus		
	7500, 7500 Fast	qMax Green Low Rox	qMax Probe Low Rox
	ViiA7™		
	QuantStudio™ 12K Flex, QuantStudio 3		
BioRad®	CFX96™	qMax Green No Rox	qMax Probe No Rox
	CFX384™		
	iCycler™, MyIQ™	qMax Green No Rox	qMax Probe No Rox
	IQ5™		
	Opticon™, Opticon™ 2	qMax Green No Rox	qMax Probe No Rox
	MiniOpticon™, Chromo4™		
Roche®	LightCycler™ 96	qMax Green No Rox	qMax Probe No Rox
	LightCycler™ 480		
	LightCycler™ Nano		
Qiagen® (Corbett)	Rotor-Gene™ Q	qMax Green No Rox	qMax Probe No Rox
	Rotor-Gene™ 3000		
	Rotor-Gene™ 6000		
Eppendorf®	Mastercycler™ ep realplex	qMax Green No Rox	qMax Probe No Rox
	Mastercycler™ ep realplex 2S		
Agilent® (Stratagene)	MX4000™	qMax Green Low Rox	qMax Probe Low Rox
	MX3000P™		
	MX3005P™		
	AriaMx		
Cepheid®	SmartCycler™	qMax Green No Rox	qMax Probe No Rox
Fluidigm®	BioMark™	qMax Green Low Rox	qMax Probe Low Rox
Illumina	Eco qPCR		qMax Probe No Rox