

# Description of Dilator Slides


## DILATOR STORAGE & DRYING CART for Savary-Galliard® Dilators

Each Dilator Cart includes four slides with membrane holders which secure Savary-Gilliard® Dilators of varying sizes, with a capacity of six dilators per slide.

The “A” slide holds six dilators between the diameter of 5-7 mm (15-21 French).

The “B” slide holds six dilators between the diameter of 8-13 mm (24-39 French).

The two “C” slides each hold six dilators (twelve dilators total) between the diameter of 14-20 mm (42-60 French).



“A” Slide  
5-7 mm  
15-21 Fr

“B” Slide  
8-13 mm  
24-39 Fr

“C” Slide  
14-20 mm  
42-60 Fr

“C” Slide  
14-20 mm  
42-60 Fr



This is the standard configuration for all dilator carts. If you require a custom configuration, please contact us.