

DANi™ Adult Patient Simulator

Product Name & Overview	DANi Adult Advanced Training Solution
Addresses Training Needs/ Curriculum	<ul style="list-style-type: none"> • Basic Life Support (BLS) • Advanced Cardiac Life Support (ACLS)
Configuration	Configurable to fit unique training needs and budgets A range of accessories available to maximize versatility
Consumables	<p>IV arm:</p> <ul style="list-style-type: none"> • Skins • Veins kit • IM pads • Fluid supply bag • Simulated blood <p>Urinary catheterization:</p> <ul style="list-style-type: none"> • Male genitalia foreskins • Foley catheter • Urinary catheterization bladders • Simulated urine <p>Additional consumables:</p> <ul style="list-style-type: none"> • Lung bags • IM pads for shoulder, thighs and buttocks • Stomach/stoma/gastronomy bags • Simulated blood • Lubricant • Stomach bags • Enema reservoir
Product Codes/SKUs	<p>Patient Care Trainer with Chest Compression</p> <ul style="list-style-type: none"> • LF04550 (Light skin tone) • LF04550M (Medium skin tone) <p>Patient Skills Trainer</p> <ul style="list-style-type: none"> • LF04502 (Light skin tone) • LF04502M (Medium skin tone) <p>Auscultation Patient Skills Trainer</p> <ul style="list-style-type: none"> • LF04503 (Light skin tone) • LF04503M (Medium skin tone)
Dimensions/ Weight	Dimensions: 22.25" H x 44.5" L x 25.5" W Ship Weight: 80 lbs.

Meet the New DANi™:
Realistic. Versatile.
Adult Patient Training Redefined.

Be
READY



Be **READY** for the future of clinical training.

DANi™ has evolved — redesigned to reflect the real-world needs of today's healthcare educators. Now a full-body adult patient, DANi delivers advanced training capabilities with a clean, natural look that's ready.

DANi features superior range of motion, with weight distributed to accurately represent a real patient for lifting and carrying techniques. A realistic, hands-on learning experience that is ideally suited for classroom, in-situ, or hospital instruction

Each version is affordable and durable with upgradable items.

Patient Care Trainer with Chest Compression

General patient care, patient positioning, bandaging, wound dressing, oxygen administration and compressible chest.

Patient Skills Trainer

Respiratory support, injection/IV with IM pads, intubation, CPR with quality performance metrics, non-invasive blood pressure, enema administration (female), catheterization, vaginal, pelvic and prostate examinations.

Auscultation Patient Skills Trainer

All the features of Patient Skills Trainer with added auscultation. Lung sounds at five anterior, ten posterior, and two midaxillary locations. Six anterior heart sites.

This manikin is intended to be used as an instructional aide for teaching purposes. Actual product may vary slightly from photo. Nasco Healthcare reserves the right to change product color, materials, supplies, or function as needed.



