

# CLEANING & CARE GUIDELINES

INSTRUCTIONS OF USE

## USE

---

Phillips Safety Products must be used only for their designated purposes.

## WARRANTY

---

Phillips Safety Products have a one-year warranty against manufacturer defects. Attention: This does NOT include normal wear and overloading of use.

## CLEANING

---

Use a mild detergent or cleaner to remove any dirt or debris.  
Do NOT use abrasives.

## DISINFECTING

---

It can be disinfected with most standard hospital-grade disinfectants containing QATs, IPA, or Hydrogen Peroxide.