

CLEANING & CARE GUIDELINES

INSTRUCTIONS OF USE

USE

Phillips Safety Products must be use for it the designated propous only.

WARRANTY

Phillips Safety Products have one year warranty against manufacturer defects. Attention: This does NOT include normal wear and overloading of use.

CLEANING

Use a mild detergent or cleaner to remove any dirt or debris.
Do NOT use abrasives.

DISINFECTING

It can be disinfected with most standard hospital grade disinfectants containing: QATs, IPA, or Hydrogen Peroxide.