



Magnifiers

UV Light

The ideal light for numerous applications.

Designed for a variety of uses, Burton UV (Woods) Lights are offered in models with magnifiers and white lights to best suit your needs. With handheld flexibility, large magnification and a five-year limited warranty, this UV light provides great value for your facility.

Burton's UV lights are depended on worldwide for performance and durability.




The Right Light

Features

- Special fluorescent bulbs produce UVA at approximately 365 nm wavelength
- Available with magnifier and white bulbs
- 4-light magnifier models include optically ground and polished glass 5 diopter magnifier lens with 2.25X total power magnification

Summary of technical data	
Rated life of bulb:	6000 hours
Diopter of magnifying lens:	5
Focal length:	8" (20 cm)
Type of bulb:	UV or white
Power requirements:	120V model standard; 230/240V models available
Product warranty:	5-year limited warranty
Origin:	Assembled in USA

Note: The optical data are nominal values based on measurements calculated according to UL/IEC performance standards at 120V



Accessories	
#1007090	UV Light Wall Bracket
#0001127PK	UV Replacement Bulbs, White, 4/Box
#1003073	UV Replacement Bulbs, UV, 4/Box