



Magnifier

# UV Light

The ideal light for numerous applications.

Designed for a variety of uses, the Burton UV (Woods) Light is also available in a model with a magnifier and a model with a magnifier and white and UV bulbs. With handheld flexibility, large magnification and a five-year limited warranty, this UV light provides great value for your facility.


## Features

- Special fluorescent bulbs produce UVA at approximately 365 nm wavelength
- 4-light magnifier models include optically ground and polished glass 5 diopter magnifier lens with 2.25X total power magnification



burton®

## UV Light product dimensions



Burton's UV lights are depended on worldwide for performance and durability.

## Summary of technical data

### UV LIGHT

Optics:	5 Diopter
Viewing Area Diameter:	8" (20 cm)
Light Source:	UV or white
Rated Life of Bulb:	6,000-hour
Power:	120v
Product Warranty:	5-year limited warranty

## Ordering Information

### UV LIGHT

Model	Description
UV501	UV Light, 2 UV Bulbs
UV502	UV Light with Magnifier, 2 UV Bulbs
UV503	UV Light with Magnifier, 2 White and 2 UV Bulbs

### ACCESSORIES

Model	Description
1007090	UV Light Wall Bracket
0001127PK	UV Replacement Bulbs, White, 4/Box
1003073PK	UV Replacement Bulbs, UV, 4/Box



© 2021 Burton Medical. All rights reserved. Burton Medical reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

burton®