

XenaLux™

Comfortable, lightweight, high quality headlight that helps you go the distance.

Features:

- 225,000 lux at 16" (40.64 cm) working distance
- 300 Watt Illuminator
- Xenon lamp provides brilliant “daylight” natural light quality
- Variable light output control
- Lightweight Head Lamp

XENALUX



burton®

Summary of technical data

XENALUX™

Illuminance:	225,000 Lux at 16" (40.64 cm)
Color Temperature:	5,600 Kelvin
CRI (Color Rendering Index):	98
Light Source:	Xenon 300W Lamp
Rated Lamp Life:	1,000 hours
Illuminator Input Voltage:	100-120VAC 50/60 Hz, 220-240VAC 50 Hz
Weight Headband, Module and Cable Spot Diameter:	6.5 ounces, 184.27 grams 20mm – 110mm @ 40cm .787" — 4.331" @ 16"
Certification/Approvals:	IEC 60601-1 Edition 3.1 IEC 60601-1-2 4th Edition FDA Class I
Product Warranty:	3-year Limited Headlamp and Module 3-year Limited Illuminator 1-year 1-year Limited Fiber Optic Cable

Headlamp:

- Lightweight and suitable for daily use
- Headband with Module Weight 4.5 ounces
- Headband, Module and Fiber Optic Cable 6.5 ounce
- Module Variable spot at 16" .78" – 4.33" (20mm-110mm)

Xenon Lamp Illuminator:

- High intensity White Light
- Ergonomically designed turret includes ACMI, Olympus, Storz & Wolf ports
- Digital Display for light level intensity and lamp hour identification
- Long lamp life, typical 1,000 hour use

Application Use:

- Surgery
- ENT
- Plastic Surgery

Ordering Information

XENALUX™

Model	Description
XNI0C	XenaLux Complete Headlight System
XNI0H	XenaLux Headlight and Cable 9' Storz Connector
XNI0I	XenaLux 300 Watt Xenon Illuminator

ACCESSORIES


Model	Description
0880300	Mobile Floor Stand for Illuminator
4000446	Joysticks Non-Sterile Autoclavable, Pack of 5
4000443	300 XenaLux Replacement Lamp
4000447	XenaLux Light Source Headlight Replacement Cable



Optional mobile Floor Stand offers vertical support



Xenon Replacement Lamp Module for Xenalux



© 2024 Burton Medical. All rights reserved. Burton Medical reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

burton®