



Major Surgery Lights

# XenaLux<sup>TM</sup>

Advanced design ensures top performance and a noticeably brighter light beam.

Powered by a 300-watt Xenon source and with light delivery via multi-fiber cables, Burton's XenaLux<sup>TM</sup> OR headlamp system offers superior intensity as well as remarkable reliability and value. The headlamp's optical system has a full range of articulation with a variable spot range, removable joystick adjustment and a new and improved vertically adjustable linkage arm for added reach and greater maneuverability. Unique design and optics allow the XenaLux<sup>TM</sup> Illuminator to provide superior features.

## Headlamp Features

- Extremely bright light at working distance
- Variable spot with flawless circular beam
- Very light weight right angle headlamp
- Preferred for comfort
- Compact optical system
- Fully-adjustable, uniquely-designed linkage arm allows for vertical adjustment of the optical module
- Adapts to any Wolf port
- Attractive standard 3.5 bundle cables

## Illuminator Features

- Produces highest intensity white light of any 300-watt Xenon source
- Quiet, cool venting system
- Rotary Turret for Wolf, Storz, Olympus, ACMI and Pilling connectors
- Hours-in-use meter
- Optional mobile floorstand offers vertical adjustment on upright assembly; the illuminator's shelf includes a handle for easy positioning




## Summary of technical data

Illuminance:	120,000+ lux (11,148+ fc)
Bulb life:	500 hours
Power (light bulbs):	300-watt
Type of bulb:	Xenon
Product warranty:	3-year limited warranty
Origin:	Assembled in USA

*Note: The optical data are nominal values based on measurements calculated according to UL/IEC performance standards at 120V*



Sample Model No: XN10C = Complete XenaLux™ Headlight System



### Accessories

#0880300

Mobile Floorstand for Illuminator