



350° swivel head display



Great for doctor offices and clinics with minimal floor space, the 599KL has a high capacity with waist-high display that can swivel 350°. For providers wishing to measure height, an optional height rod is available.

Advantages

- EMR Capable**
Scale interfaces with a computer, EMR software or other electronic device via Health o meter[®] Professional Connectivity Solutions
- Value-Added Features**
Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost
- 350° Swivel Display**
Eye-level rotating display can be read from almost any angle and allows for ease of use and patient privacy
- Calculates Body Mass Index (BMI)**
Scale functions include the ability to calculate the patient's BMI



Specifications

- Capacity: 600 lb / 272 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size (w x d x h):
14 1/4" x 14 1/4" x 2 5/8"
362 mm x 362 mm x 67 mm
- Product Footprint (w x d x h):
14 3/4" x 23 1/4" x 42 1/2"
375 mm x 591 mm x 1080 mm
- Product Weight: 20 lb / 9 kg
- Height Rod: Optional 499KLROD
- Display: 1" / 25 mm LCD, 350° Swivel
- Wheels: 2
- Connectivity: USB, optional Pelstar® wireless technology
- Power Source: 100-240V adapter (ADPT50) included,
6 C-cell batteries (not included)
- Functions: LB / KG Conversion, LB / KG Lock Out , BMI, Zero, Tare,
Hold / Release, Auto Zero, Auto Off, Parent / Child Weighing Function
- 2 year limited warranty

Optional height rod
for 599KL or 599KG
(See chart below)

Model	Description
599KG	Heavy duty waist high digital scale, kilograms only
599KLHR	599KL with height rod included
599KGHR	599KG with height rod included
SS-599KL	ScaleSurance 2 year extended warranty for 599 models*
499KLROD	Optional height rod for 599KL & 599KG
ADPT50	Adapter included, item # for replacement or spare
PLUGSET50	Plug adapters for use outside North America
PELSTARONE-P	Wireless thermal printer

* Not available in all countries