

LID PLATE

PRODUCT INFORMATION

PRODUCT INFORMATION

LID PLATE | LS-JLP-BK, LS-JLP, LS-KJLP-S, LS-KJLP-L



The Lid Plate is a cutting-edge instrument designed to deliver unparalleled stability and safeguard the delicate ocular globe during lid dissection procedures. For this purpose, the instrument features a wide, concave blade that follows the anatomy of the globe. Crafted with precision, this indispensable tool is available in four models to accommodate diverse surgical preferences. These include the Plain Stainless, Autoclavable Jaeger Lid Plate, Khan-Jaeger Lid Plate, and Jaeger Lid Plate, each catering to diverse surgical preferences and ensuring optimal performance in the operating room.

PRODUCT SPECIFICATION

MATERIAL: Stainless Steel

COLOR: Gray, Black