

PTHC65GUC



34.25" x 23.5" x 23.5" (H x W x D)

Undercounter single chamber passive heating blanket warmer with glass door, stainless steel interior and exterior, digital thermostat, and lock

Highlights:

Designed and assembled in the USA

Temperature range: 90° to 140°F

No internal fans or moving parts ensure silent, safe, and efficient operation

Product Features:

| | |
|---------------------------------|--|
| Undercounter Warming Cabinet | 24" wide blanket warmer designed and constructed in the USA |
| 90 to 140°F temperature range | Designed and set to keep blankets warm in hospitals, medical and physical therapy offices, spas, and other facilities |
| Passive heating | Fixed heaters located at the bottom of the cabinet with no fans or moving parts to minimize lint buildup for safer, more efficient use |
| Built-In capable | Front-breathing design lets you make the best use of space by installing your appliance under the counter |
| Digital thermostat | Electronic controls are conveniently located externally for easy temperature adjustment |
| Electronic readout | Large LED temperature display ensures your staff can easily monitor the interior temperature |
| Silent operation | No moving parts to ensure a quiet performance |
| Double pane tempered glass door | Heat-safe glass door offers a full view of stored contents, with a stainless steel door trim to highlight the interior |
| Professional handle | Stainless steel handle completes the professional look |

PTHC65GUC Specifications:

| Overview | |
|----------------------------|---------------------|
| Nominal Width | 24.0" (61 cm) |
| Nominal Depth | 24.0" (61 cm) |
| Nominal Height | 35.0" (89 cm) |
| Height of Cabinet | 34.25" (87 cm) |
| Width | 23.5" (60 cm) |
| Depth | 23.5" (60 cm) |
| Depth with Handle | 25.38" (64 cm) |
| Depth with door at 90° | 45.75" (116 cm) |
| Capacity | 4.57 cu.ft. (129 L) |
| Door | Glass |
| Cabinet | Stainless Steel |
| US Electrical Safety | ETL |
| Canadian Electrical Safety | ETL-C |
| Amps | 5.5 |
| Voltage/Frequency | 115 V AC/60 Hz |
| Weight | 125.0 lbs. (57 kg) |
| Parts & Labor Warranty | 1 Year |
| UPC | 761101108056 |
| Warming Cabinet | |
| Door Swing | RHD |
| Reversible | Factory Reversible |
| Adjustable Shelves | Yes |
| Shelf Type | Chrome |
| Shelf Qty | 3 |
| Thermostat Type | Digital |
| Interior Height | 26.0" (66 cm) |
| Interior Width | 19.38" (49 cm) |
| Interior Depth | 17.38" (44 cm) |
| Level Legs | 4 |
| Dolly | Optional |
| Temperature Range | 90 to 140°F |

