



PARAGON

DISPOSABLE SINGLE USE HOSE COVER

To install

1. Take the end of the hose with male connectors and insert it into the hose cover. Gently maneuver the hose through until the hose is completely covered.
2. Connect the male connectors from the hose into the female connectors on the cooling system. An audible click can be heard when the connectors are properly connected.

To remove

1. Hold the end of the hose with female connectors and pull the hose cover off from the opposite end of the hose.
2. Dispose of the used hose cover.



PARAGON

DISPOSABLE SINGLE USE HOSE COVER

To install

1. Take the end of the hose with male connectors and insert it into the hose cover. Gently maneuver the hose through until the hose is completely covered.
2. Connect the male connectors from the hose into the female connectors on the cooling system. An audible click can be heard when the connectors are properly connected.

To remove

1. Hold the end of the hose with female connectors and pull the hose cover off from the opposite end of the hose.
2. Dispose of the used hose cover.