

PTHC125G



60" x 23.63" x 23.75" (H x W x D)

Single chamber passive heating blanket warmer with a self-closing glass door, stainless steel interior, digital thermostat, stainless steel cabinet, and lock

Highlights:

Designed and assembled in the USA

No internal fans or moving parts ensure silent, safe, and efficient operation



Stainless steel construction inside and out for durable use and easy cleaning

Product Features:

90 to 140°F temperature range	Designed and set to keep blankets warm in hospitals, medical and physical therapy offices, spas, and other facilities
Passive heating	Fixed heaters located at the bottom of the cabinet with no fans or moving parts to minimize lint buildup for safer, more efficient use
Digital thermostat	Electronic controls are conveniently located for easy temperature adjustment
Electronic readout	Large LED temperature display ensures your staff can easily monitor the interior temperature
Silent operation	No moving parts to ensure a quiet performance
Double pane tempered glass door	Heat-safe glass door offers a full view of stored contents, with a stainless steel door trim to highlight the interior
Professional handle	Stainless steel handle completes the professional look
Self-closing door	Door will automatically close when left ajar, preventing items from accidentally being left exposed to ambient temperatures
Stainless steel interior	Interior walls are constructed from stainless steel to ensure lasting durability and easy cleaning

PTHC125G Specifications:

Overview	
Nominal Width	24.0" (61 cm)
Nominal Depth	24.0" (61 cm)
Nominal Height	60.0" (152 cm)
Height of Cabinet	60.0" (152 cm)
Width	23.63" (60 cm)
Depth	23.75" (60 cm)
Depth with Handle	25.63" (65 cm)
Depth with door at 90°	46.0" (117 cm)
Capacity	10.57 cu.ft. (299 L)
Door	Glass
Cabinet	Stainless Steel
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	5.5
Voltage/Frequency	115 V AC/60 Hz
Weight	160.0 lbs. (73 kg)
Parts & Labor Warranty	1 Year
UPC	761101066127
Warming Cabinet	
Door Swing	RHD
Reversible	No
Adjustable Shelves	Yes
Shelf Type	Chrome
Shelf Qty	4
Thermostat Type	Digital
Interior Height	45.25" (115 cm)
Interior Width	20.0" (51 cm)
Interior Depth	16.5" (42 cm)
Dolly	Optional
Temperature Range	90 to 140°F



Puretherm Warming Cabinet