

THE 245EHR-3105

DIGITAL HEIGHT ROD FOR THE 3105 SERIES OF SCALES

The Health o meter® Professional digital 245EHR-3105 is an optional state-of-the-art digital height rod for the Antimicrobial 3105 series of scales. The market-leading height rod provides the most accurate height measurement in its class, with a graduation/resolution of 1 mm/0.1" and a full measuring range of 11 ¾" – 79 ½" / 30cm – 202cm and a digital range of 3ft 11.2" – 6ft 7.5" / 120 cm - 202 cm.



Built-In Reset Feature

Proprietary technology eliminates mechanical sensors, reducing failures, increasing accuracy and reliability

The new accessory uses solid-state technology, which provides a facility confidence that a patient's height will be measured accurately, and workflow will not be interrupted. The built-in reset feature increases accuracy and reliability, and the folding headpiece limits the height rod's profile, providing safety to patients and caregivers. The optional digital height rod has quick and easy to assemble brackets.



TECHNICAL SPECIFICATIONS

Digital Height Range

3ft 11.2" – 6ft 7.5" / 120 cm - 202 cm
Resolution: 1 mm / 0.1"

Mechanical Height Range

11 ¾" – 79 ½" / 30cm – 202cm
Graduation: 1 mm / 1/16 inch

Product Footprint with 3105KL-AM

(height rod extended)

23 ¼" x 28" x 80"
(591 mm x 711 mm x 2032mm)

Product Footprint with 3105KL-AM

(height rod in home position)

23 ¼" x 28" x 52 ¼"
(591 mm x 711 mm x 1327mm)

Product Weight with 3105KL-AM

90 lb / 41 kg

Warranty

2-year limited manufacturer warranty

245EHR-3105 Compatible Scales

3105KL-AM
3105KG-AM
3105KL-AM-BT
3105KG-AM-BT

ORDERING INFORMATION

MODEL	DESCRIPTION
245EHR-3105	Optional digital height rod for 3105 Series
SS-245EHR	ScaleSurance® 2 year extended warranty*

* Not available in all countries

Health o meter®
Professional

Health o meter® Professional 245EHR-3105

weigh**easier**®