

## THE 3001 HR

# MECHANICAL HEIGHT ROD FOR THE 3001 SERIES OF SCALES

The Health o meter® Professional 3001HR is an optional mechanical height rod for the Antimicrobial 3001 series of scales. The durable aluminum height rod withstands frequent use, and rigorous laboratory testing shows that the 3001HR lasts 2 ½ times longer than similar products.

The easy-to-use smooth sliding headpiece and telescoping rod allow the caregiver to read the height measurement quickly and easily. The folding headpiece limits the height rod's profile, providing safety to patients and caregivers. The 3001HR has an extensive measuring range of 2"-90" in 1/16" increments (6 cm-230cm in 1mm increments). The optional mechanical height rod has upgraded brackets making assembly quick and easy.



## TECHNICAL SPECIFICATIONS

### Measuring Range

2" – 90" / 6cm – 230cm

Graduation: 1/16" / 1mm

### Warranty

2-year limited manufacturer warranty

### 3001HR Compatible Scales

3001KL-AM

3001KG-AM

3001KL-AM-BT

3001KG-AM-BT

## ORDERING INFORMATION

MODEL	DESCRIPTION
3001HR	Optional mechanical height rod for 3001 Series

**Health o meter®**  
P r o f e s s i o n a l

Health o meter® Professional 3001HR

weigh**easier**®