

<b>Product Name:</b>	<b>Examination Gloves</b>
<b>Product Application:</b>	<b>NEST gloves can be used in household, industrial, medical, beauty and other industries, and it is an essential hand protection product.</b>
<b>Product Materials:</b>	<b>Latex / Nitrile + colloidal oatmeal</b>
<b>product Specifications:</b>	Temperature range: stored at room temperature  Shelf life: 3 Years after date of production ( ensure package is in good )
<b>Sterilization:</b>	No
<b>BSE/TSE:</b>	These products are deemed animal free.
<b>Performance Testing:</b>	Each manufacturing lot is sampled and tested in accordance with standard operating procedures.  Appearance inspection: qualified                      Tensile Strength Test: qualified Leak test: qualified    Packaging inspection: qualified
<b>Features:</b>	<ol style="list-style-type: none"> <li>1. Single use only</li> <li>2. Made in Malaysia</li> <li>3. Ambidextrous: Can be worn left and right, regardless of left or right hand</li> <li>4. Glove mouth rolled edge design, easy to wear, prevent rollback</li> <li>5. Textured outer surface of the glove increases friction and facilitates gripping</li> <li>6. Powder-Free</li> <li>7. Non-Sterile</li> <li>8. Boxed for easy extraction.</li> <li>9. <b>Available in three materials</b> <ul style="list-style-type: none"> <li><b>(1) Nitrile Examination Glove:</b> <ul style="list-style-type: none"> <li>• Tensile Strength: Minimum 18.0 MPa (Before Aging)</li> <li>• Palm thickness: 0.07 mm</li> </ul> </li> <li><b>(2) Latex Examination Glove:</b> <ul style="list-style-type: none"> <li>• Natural rubber latex and compounding materials</li> <li>• Tensile Strength: Minimum 18.0 MPa (Before Aging)</li> <li>• Palm thickness: 0.11 mm</li> </ul> </li> <li><b>(3) Colloidal Oatmeal Coated Nitrile Glove</b> <b>Colloidal Oatmeal Coated Latex Glove</b> <ul style="list-style-type: none"> <li>• Oats Extractions</li> </ul> <p>Oats extractions is a patented coating recognized by the FDA as a skin protectant, which have been proven by clinical studies could reduce and control rash and dry skin symptoms. Clinical studies have shown that colloidal oatmeal could reduce and control rash, dry skin.</p> </li> </ul> </li> </ol>

**Product Range**

**Nitrile Examination Glove**

Cat. No.	Size	Palm Width (mm) ±4	Length	/Inner Box	/Case
902001	XS	76	Min 230mm (ASTM)	100	1000
902011	S	86		100	1000
902021	M	98		100	1000
902031	L	107		100	1000
902041	XL	115		90	900

**Latex Examination Glove**

Cat. No.	Size	Palm Width (mm) ±4	Length	/Inner Box	/Case
901001	XS	76	Min 230mm (ASTM)	100	1000
901011	S	86		100	1000
901021	M	98		100	1000
901031	L	107		100	1000
901041	XL	115		90	900

**Colloidal Oatmeal Coated Nitrile Glove**

Cat. No.	Size	Palm Width (mm) ±4	Length	/Inner Box	/Case
903011	S	84	Min 230mm (ASTM)	100	1000
903021	M	94		100	1000
903031	L	105		100	1000

**Colloidal Oatmeal Coated Latex Glove**

Cat. No.	Size	Palm Width (mm) ±4	Length	/Inner Box	/Case
904011	S	84	Min 230mm (ASTM)	100	1000
904021	M	94		100	1000
904031	L	105		100	1000

Product Picture

Nitrile Examination Glove



Due to chromatic aberration of the display screen, there may be chromatic aberration between the picture and the real object.

Latex Examination Glove



The picture is only for reference, please make the object as the standard.