

**TABLE TOP CASSETTE HOLDER
TT-CSTHDR**

PRODUCT INFORMATION

COPYRIGHT ©2025. ALL RIGHTS RESERVED



**PHILLIPS SAFETY
PRODUCTS INC.**

PRODUCT INFORMATION
TABLE TOP CASSETTE HOLDER | TT-CSTHDR



The Table Top Cassette Holder is designed for lateral imaging. Crafted from top-grade steel, this holder guarantees enduring performance and durability. Plus, it is engineered for stability, featuring a 15" x 15" base that effortlessly slides beneath patients and pads. Featuring a 1.25" channel depth, this Table Top Cassette Holder accommodates a variety of popular imaging CR/DR Cassettes and Grid Cassettes. Adjustments are seamless, facilitated by the smooth-sliding clamps on the 24" steel pole. In addition, the receptor handles snugly fit into the upper clamp, securely cradled on the expansive base. Each clamp is lined with anti-slip soft rubber for gentle handling and protection of the receptor at all times.

SPECIFICATION

COLOR: Gunmetal

OVERALL SIZE: 15"L x 15"W x 27"H

FEATURES: Adjustable Height

CHANNEL DEPTH: 1.25"

ADJUSTMENT HEIGHT RANGE: 1" to 23.5"