

# THE 4060

## ONE-WAY SLIDE CUSHION REPOSITIONING AID

The One-Way Slide Cushion was designed to safely position patients/ residents in chairs, wheelchairs, beds, or cars reducing the risks of slips and falls. The One-Way Slide Cushion uses a low-friction tubular design that only slides to the back of the chair, not forward. This prevents the patient/resident from egressing forward out of the chair. The cushion simplifies caregiver assistance, eliminating the need for lifting, while patients can easily adjust it independently.

### PRODUCT SPECIFICATIONS

**Product Size**  
18" x 16"

**Material**  
Polyurethane coated polyester, Latex free

**Weight Limit**  
550 lb/ 250 kg

**Machine-Washable**  
Can easily be washed under water temperature 160° F/71° C. Do not tumble dry, do not iron, do not bleach



- Low friction cushion that holds patient/resident from sliding forward in their chair
- Easily adjust patient/resident back in their seat to maintain posture
- Patient/ resident maintains a safe, secure posture
- Prevents fall-related injury
- Clearly labeled for proper orientation
- Soft, comfortable, breathable/ wicking material; tubular design
- Easy for patients/residents to do by themselves
- Easy for staff to help without having to lift
- Can be used in a car, chair, bed, or underneath any pressure relief cushion or incontinent pad

This device prevents the patient/ resident from sliding forward and falling.



### ORDERING INFORMATION

MODEL	DESCRIPTION	
4060	One-Way Slide Cushion Repositioning Aid	1/ea