





fitted with	1 x LED module
luminous power	50'000 Lux / 0.5 m
Illuminancefield 0.5 m color rendering index (CRI)	Ø 18 cm (7.1 in) > 93
color temperature work	4400 K
equipment	electronic converter
connected load	100-240 V; 50/60 Hz
mains lead	1.8 m (5.9 ft) hospital grade plug
luminaire body	PMMA
weight (net)	1.1 kg (2.4 lbs)
power consumption	13 W
usage	rocker switch (I/O) at the arm
class of protection	II
accessory	wall bracket
norms	EN 60601-1, EN 60601-2-41
special features	positioning handle
certifications	ETL/cETL
system of protection	IP 20
order number	D15757120 - pure white



light distribution







fitted with	1 x LED module
luminous power	50'000 Lux / 0.5 m
Illuminancefield 0.5 m color rendering index (CRI)	Ø 18 cm (7.1 in) > 93
color temperature work	4400 K
equipment	electronic converter
connected load	100-240 V; 50/60 Hz
mains lead	1.8 m (5.9 ft) hospital grade plug
luminaire body	PMMA
weight (net)	1.1 kg (2.4 lbs)
power consumption	13 W
usage	rocker switch (I/O) at the arm
class of protection	II
accessory	rail clamp
norms	EN 60601-1, EN 60601-2-41
special features	positioning handle
certifications	ETL/cETL
system of protection	IP 20
order number	D15757150 - pure white



light distribution







fitted with	1 x LED module
luminous power	50'000 Lux / 0.5 m
Illuminancefield 0.5 m color rendering index (CRI)	Ø 18 cm (7.1 in) > 93
color temperature work equipment	4400 K electronic converter
connected load	100-240 V; 50/60 Hz
mains lead	1.8 m (5.9 ft) hospital grade plug
luminaire body	PMMA
weight (net)	1.1 kg (2.4 lbs)
power consumption	13 W
usage	rocker switch (I/O) at the arm
class of protection	II
accessory	table clamp
norms	EN 60601-1, EN 60601-2-41
special features	positioning handle
certifications	ETL/cETL
system of protection	IP 20
order number	D15757100 - pure white



light distribution





fitted with	1 x LED module
luminous power	50'000 Lux / 0.5 m
Illuminancefield 0.5 m color	Ø 18 cm (7.1 in)
rendering index (CRI)	> 93
color temperature work	4400 K
equipment	electronic converter
connected load	100-240 V; 50/60 Hz
mains lead	1.8 m (5.9 ft) hospital grade plug
luminaire body	PMMA
weight (net)	1.1 kg (2.4 lbs)
power consumption	13 W
usage	rocker switch (I/O) at the arm
class of protection	II
accessory	roll stand



norms	EN 60601-1, EN 60601-2-41
special features	positioning handle
certifications	ETL/cETL
system of protection	IP 20
order number	D15757110 - pure white

light distribution

