



### Advantages

- Easy to Move**  
Lightweight frame with 4 casters make portability easy for bringing the scale to the patient's location.
- High Capacity**  
High 600 lb / 270 kg capacity provides flexibility to weigh almost any patient.
- No Special Shipping Requirements**  
The 594KL can be shipped via ground service, no need for a pallet or loading dock
- Designed for Patient Comfort**  
Large seat with moving armrests allow for easy access on and off the scale.
- Convenient Storage Features**  
Folding seat back and armrests make the scale simple to store.
- Value-Added Features**  
Operating functions (LB/KG Conversion, LB/KG Lock Out, Body Surface Area (BSA), BMI, Zero, Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost



**Specifications**

- Capacity: 600 lb / 270 kg
- Resolution: 0.2 lb / 0.1 kg
- Seat Size (w x d):  
18" x 18 ½"  
457 mm x 470 mm
- Product Footprint (w x d x h):  
23 ¼" x 40 ⅝" x 37 ⅞"  
590 mm x 1020 mm x 960 mm
- Product Weight: 39 lb / 18 kg
- Display: 1 ⅜" / 35 mm LCD
- Wheels: 4 Casters for ease of movement
- Power Source: 100-240V adapter (69-00001) included, 6 AA batteries not included
- Functions: LB/KG Conversion, LB/KG Lock Out, BSA, BMI, Zero, Tare, Hold/Release, Auto Zero, Auto Off,
- 2 year limited warranty

Model	Description
SS-594KL	ScaleSurance 2 year extended warranty*
69-00001	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America