



fitted with	LED - Energy efficiency category A+
color rendering index (CRI)	> 95
color temperature	6500 K
work equipment	electronic converter
mains lead	59 in (1.5 m)
lens; dimensions	Ø 6.3 in (160 mm)
lens; diopters	1.4 in (3.5 cm)
luminaire body	polycarbonate & acrylonitrile butadiene styrene copolymer (PC & ABS)
weight (net)	1.3 lbs (0.6 kg)
power consumption	7 W
usage	wireless switch
class of protection	II
norms	EN 60601-1
special features	battery, ergonomic handle, scratch-proof plastic lens, protective cover, continuously dimmable, switchable wood-light
order number	D16067000 - pure white



light distribution

