

**Health o meter**<sup>®</sup>  
P r o f e s s i o n a l



### Advantages

- ✓ **Measuring Tape Included**  
Easy-to-read measuring tape allows user to measure both weight and height simultaneously
- ✓ **Dual Unit Measurement Bar**  
The mechanical poise bars allow for easy reading of weight in lb or kg

**Health o meter**<sup>®</sup>  
P r o f e s s i o n a l  
www.homscales.com

# Health o meter<sup>®</sup>

P r o f e s s i o n a l



Easy to read  
poise bars

### Specifications

- Capacity: 135 lb / 62 kg
- Graduation: 1 oz / 20 g
- Tray Size: 24" (w) x 16 1/2" (d) x 2 7/8" (h), 610 mm x (w) x 420 mm (d) x 70 mm (h)
- Easy to read poise bars for dual measurement in lb and kg
- Easy to clean plastic tray
- Measuring Tape: 0" - 23 3/4" / 0 cm - 60 cm
- 2 Year Limited Warranty

**Health o meter<sup>®</sup>**  
Professional  
[www.homscales.com](http://www.homscales.com)