

**Examination Luminaire**  
**HALUX N30-1 F**  
 # D16064000



<b>fitted with</b>	LED module - energy efficiency category A+
<b>illuminance field @ 50cm</b>	Ø 18cm
<b>illuminance @ 50cm</b>	30,000 lx
<b>color rendering index (CRI)</b>	> 93
<b>color temperature</b>	4400K
<b>work equipment</b>	electronic plug-in transformer
<b>connected load</b>	100-240 V; 50/60 Hz
<b>mains lead</b>	approx. 6.6 ft (2 m)
<b>luminaire body</b>	polycarbonate & acrylonitrile butadiene styrene copolymer (PC & ABS)
<b>weight (net)</b>	approx. 9.7 lbs (4.4 kg)
<b>power consumption</b>	8W
<b>usage</b>	keypad at luminaire head
<b>balance of articulation joints</b>	friction
<b>class of protection</b>	II
<b>norms</b>	EN 60601-1; EN 60601-2-41
<b>special features</b>	continuously dimmable, maintenance free
<b>fastening</b>	rolling stand
<b>order number</b>	D16064000 - pure white

light distribution

