

# NEW BowTie™ Bags are Even Bigger!

## Drawstring Bags for Large Biohazard XL Bins

- Instant closure - no wire tie needed
- NEW 30 gallon or 12 gallon capacities
- Safely & easily remove from bin or autoclave
- Color changing sterilization indicator
- Hands-free, foot pedal operated lid on bin



↑ More than double the capacity of existing bags ↑



**A8003R**  
30 Gallon Bag  
33x41in.

Enjoy the drawstring convenience of BowTie™ Biohazard Bags in a new larger, 30 gallon capacity!

The exclusive drawstring loops on BowTie™ Biohazard Bags provide instant-cinch closure on your waste materials and a handle to easily remove bags from their bins. They also allow venting during sterilization and provide an easy, safe way to remove bags from a hot autoclave.

BowTie™ Bag's 2.25 mil thick HDPE is suitable for autoclaving in 121°C cycles.

BowTie™ Bags are printed with a number of identifying markers such as large biohazard symbols and temperature limit warnings.

They are also outfitted with a large marking area and a sterilization patch that changes color after autoclaving



MTC's BowTie™ Bags fit their respective biohazard bins perfectly. The hands-free lid offers quick disposal and a cleaner, touchless experience with just a simple step on the pedal.

Whether you need a small or large capacity bag, BowTie™ Bins and Bags are a convenient addition to your hazardous waste management protocol.

**A8000B bin couples with  
12 gallon BowTie™ bags**

**A8002B XL bin is an ideal fit for  
30 gallon BowTie™ Bags.**



**A8001R  
12 Gallon**

**A8000B  
10x14x19"**



**A8003R  
30 Gallon**

**A8002B  
19x19x33"**

**BowTie™ Biohazard Bags - w/ Drawstring, Marking Area, Sterilization Indicator**

**New!**

Cat No.	Capacity	Dimensions	Bin	Qty
A8003R	30 Gallon (113L)	33x41" (84x104cm)	A8002B	100
A8001R	12 Gallon (45L)	25x35" (64x89cm)	A8000B	100

**BowTie™ Biohazard Bins - w/ Attached Lid and Hands-Free Foot Pedal**

Cat No.	Capacity	Dimensions	Bag	Qty
A8002B	26 Gallon (100L)	19x19x33" (48x48x84cm)	A8003R	1
A8000B	8 Gallon (30L)	10x14x19" (25x36x48cm)	A8001R	1



Bags are slightly larger than bin for a secure fit when folded around the rim