

# Protective Reservoir Lids

## Now *FREE* with 50mL MTC Reservoirs



- 5 Lids included *FREE* with all standard MTC 50mL & 25mL reservoirs
- Protect solutions from airborne contaminants
- Inhibit evaporation
- Designed to be easily pulled apart when stacked
- Supplied in sterile bags of 5



### Keep Your Samples Safe



Expanding on our popular 25mL reservoir lids, new 50mL lids provide security and protection for valuable reagents. These polystyrene lids perfectly line the perimeter of MTC's 25mL and 50mL reservoirs, providing complete protection. Reservoir lids ensure the security of solutions by inhibiting evaporation and protecting from airborne contaminants.

Reservoir lids are RNase, DNase, and pyrogen-free. They are packaged in sterile bags of 5.

### Simple and Convenient Design

These lids are designed with the same approach to convenience as the reservoirs they pair with. The curved corner design allows lids to cover reservoirs with a perfect fit. At the same time, these lids are built to be easily removed both from their reservoir and from their stack of lids.

MTC's lids are easy to handle and won't risk being damaged by a user struggling to remove one from the top of the stack.

