



mayo stands

foot-pedal operated

- Base designed to move under low clearance equipment with maximum stability
- Convenient foot-pedal operated lock adjusts height from 36" to 62"
- **co•re casters™** corrosion-resistant, maintenance-free
- All welded stainless steel construction



P-1065-SS

- Removable extra large 20" x 25" stainless steel tray for easy cleaning
- Two 2" corrosion-resistant Gemini **co•re casters™**
- 50 lbs. weight capacity evenly distributed



P-1066-SS

- Removable 16¼" x 21¼" stainless steel tray for easy cleaning
- Two 2" corrosion-resistant Gemini **co•re casters™**
- 50 lbs. weight capacity evenly distributed



P-1068-SS

- Removable 12⁵/₈" x 19¹/₈" stainless steel tray for easy cleaning
- Two 2" corrosion-resistant Gemini **co•re casters™**
- 50 lbs. weight capacity evenly distributed



P-1069-SS

- Removable large 16¼" x 21¼" stainless steel tray for easy cleaning
- Four 2" corrosion-resistant Gemini **co•re casters™**
- 50 lbs. weight capacity evenly distributed



P-1066-A-SS

- Removable 16 $\frac{1}{4}$ " x 21 $\frac{1}{4}$ " stainless steel tray for easy cleaning
- Adjust height 36" to 60" - friction lock
- Two 2" corrosion-resistant Gemini **co•re casters**[™]
- 50 lbs. weight capacity evenly distributed



P-1068-A-SS

- Removable 12 $\frac{5}{8}$ " x 19 $\frac{1}{8}$ " stainless steel tray for easy cleaning
- Adjust height 36" to 60" - friction lock
- Two 2" corrosion-resistant Gemini **co•re casters**[™]
- 50 lbs. weight capacity evenly distributed



mayo stands

hand operated

- Base designed to move under low clearance equipment with maximum stability
- **co•re casters**[™] corrosion-resistant, maintenance-free
- All welded stainless steel construction



P-1069-A-SS

- Removable large 16 $\frac{1}{4}$ " x 21 $\frac{1}{4}$ " stainless steel tray for easy cleaning
- Adjust height 36" to 60" - friction lock
- Four 2" corrosion-resistant Gemini **co•re casters**[™]
- 50 lbs. weight capacity evenly distributed



P-66

- Removable 19 $\frac{1}{8}$ "L x 12 $\frac{5}{8}$ "W stainless steel tray on powder-coated frame
- Adjust height from 32" to 53" - friction lock
- Two 3" corrosion-resistant **co•re casters**[™]
- 20 lbs. weight capacity evenly distributed



P-65 Medical Tray Unit

- Removable 19 $\frac{1}{8}$ "L x 12 $\frac{5}{8}$ "W stainless steel tray on chrome plated base
- Swivel height adjust from 30" to 40"
- Four 3" corrosion-resistant **co•re casters**[™]
- 40 lbs. weight capacity evenly distributed