

LEAD TABLE SHIELD VELCRO STYLE
LB-LTSV-1218 | LB-LTSV-1224 | LB-LTSV-1230 | LB-LTSV-1236

PRODUCT INFORMATION

COPYRIGHT ©2025. ALL RIGHTS RESERVED



PRODUCT INFORMATION

LEAD TABLE SHIELD VELCRO STYLE | LB-LTSV-1218 | LB-LTSV-1224 | LB-LTSV-1230 | LB-LTSV-1236



The Lead Table Shield Velcro Style are easy-to-use leaded panels that help decrease exposure to dangerous scatter radiation in any virtual imaging environment. This Leaded Table Shield features a unique modular design that ensures a perfect fit no matter the imaging table's size. Each leaded panel is 0.5mm Pb equivalent lead protection for the lower body. The Lead Table Shields features premium industrial hook & loop strips for quick installation and superior attachment to the imaging table. Available in Blue and Black colors and different sizes: 18" x 12", 24" x 12", 30" x 12", and 36" x 12".

SPECIFICATION

SIZES AVAILABLE: 18" x 12", 24" x 12",

30" x 12", and 36" x 12"

PROTECTION LEVEL: 0.50mm Pb LE

MATERIAL: Lead

COLOR: Black or blue

FEATURES: Premium industrial hook & loop strips

PRODUCT INFORMATION