

# CADDY SHIELDS

PRODUCT INFORMATION



**PHILLIPS SAFETY**  
PRODUCTS INC.

## PRODUCT INFORMATION

### CADDY SHIELDS



The Caddy Shields is a device designed to carry the weight of your apron by standing behind the device while the apron is draped over our proprietary designed hanger. This allows for you to have X-Ray protection without bearing the weight of the apron. Our H brace design allows you to stand as close to the apron as possible while allowing your arms and hands to move freely. The Caddy Shields moves effortlessly to and from your desired position. Once in position, we offer 4 locking casters. The aluminum ergonomic design makes the Caddy easy to move with minimal effort. Additionally, it features an adjustable height of 47" up to 69", making it suitable for users of varying heights.

## SPECIFICATION

---

**BRAND:** Phillips Safety

**PRODUCT TYPE:** Lead Apron Rack

**RACK TYPE:** Mobile Rack

**MATERIAL:** Aluminum

**FEATURES:** H brace design, 4 locking casters, Adjustable Height of 47" up to 69".