

**VETERINARY RADIATION HAND PROTECTION –
OPEN PALM MITT**

PRODUCT INFORMATION

PRODUCT INFORMATION

VETERINARY RADIATION HAND PROTECTION – OPEN PALM MITT | RM-OP-0.5



The Veterinary Radiation Hand Protection – Open Palm Mitt is a protective accessory designed to shield the back of the hand. It features a hook & loop strap closure across the palm for quick and easy removal. Offering 0.5mm Pb Lead Equivalency protection, the mitts are sold as a pair.

SPECIFICATION

COLOR: Black

QUANTITY: Pair

PROTECTION: 0.5mm Pb Lead Equivalency

FEATURES: Adjustable Fit, Hook and Loop Strap Closure

SIZE: 13" from top to bottom; 9" across

PRODUCT INFORMATION