

Cabinet Style, Space Saver Table with Adjustable Backrest and Leg Rest



MODEL 91002-27

- No additional charge for adjustable backrest
- All laminated construction
- Backrest has 6 settings of approximately: 60°, 55°, 50°, 43°, 36° & 27°
- Leg rest has 4 settings of approximately: 0°, 12°, 24° & 36°
- Two, closed-door, storage compartments with access from both sides
- 1 adjustable shelf behind each door
- Seamless rounded corner top
- 2" firm foam padding
- Heavyweight, knit-backed upholstery
- Ships assembled
- 350 lbs. load capacity under normal use (159 kg)

SPECIFICATIONS

91002-27	54"/72" Extended L x 27" W x 31" H
Top	54" L x 27" W
Leg Rest	18" L x 22" W


Prop 65 Warning



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1

Options



Laminate Colors



Upholstery Colors

