

Open Base, Power, Changing Table

MODEL 80110-72

- Features easy-clean, one-piece top
- Two lock-in-place, padded side rails
- Four, smooth-roll, locking, 3" casters
- Electric power height adjustment with foot control
- Ultra-strong, uni-frame base and top construction
- Heavy duty, all welded, lift mechanism
- Durable gray powder-coated steel frame
- Self-lubricating pivot points
- Smooth 110 volt electric height adjustment
- Hands-free low voltage foot control
- Easy service, self contained, electrical components
- Seamless, rounded corner, upholstered top
- 2" (5 cm) ultra-firm padding
- Four, locking casters
- Low height may help achieve ADA tax credit
- Made in the USA by Clinton Industries
- 400 lbs. (181 kg) load capacity under normal use



80110-72 SPECIFICATIONS

Overall 72" L x 27" W x 18-35" H



California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black