

Open Base, Power XL Table with One Piece Top and Casters

CLINTON OPEN BASE & SHROUDED
POWERXL
Extra width, Extra capacity, Extra Affordable

MODEL 8019987

- Features easy-clean, 30" wide, one-piece top
- Four, 3" heavy duty, rubber wheel, easy-roll, locking casters
- 600 lbs. (272 kg) load capacity under normal use
- Smooth, electric, raise and lower height adjustment
- Heavy duty, steel frame is fabricated in the USA by Clinton Industries
- Dependable scissors-lift features an H-leg design for extra stability
- Rounded corner 1" plywood top with 2" of ultra firm padding covered by heavy weight upholstery
- Harden steel bolts and self lubricating pivot bearings for years of dependable service
- Low height may help achieve ADA tax credit



8019987 SPECIFICATIONS

Overall 72" L x 30" W x 18-35" H

California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.



Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black

OPTIONS

Positioning Straps Set

- 2" x 92" strap is secured with quick hook and loop Velcro® fasteners
- 2 straps per set (92" x 2")

080



Safety Rails (pair)

- 3 height positions
- Plunger height adjustment

082



Hand-Held Controller

- Alternative to foot control
- Stores easily under top when not in use

084



030 Paper Dispenser


- Accommodates paper rolls from 18" (45.72cm) to 21" (53.34cm) wide

030

Paper Holder

- Clear plastic with elastic strap and buckle

040



Paper Dispenser/Holder Combo

- Paper dispenser and paper holder ordered as one option

034