

Open Base, Power XL Table with Adjustable Backrest, Drop Section and Casters

MODEL 8039987

- Features easy-clean, 30" wide top with adjustable pneumatic backrest and drop section
- Four, 3" heavy duty, rubber wheel, easy-roll, locking casters
- 600 lbs. (272 kg) load capacity under normal use
- Pneumatic backrest adjusts from either side and raises 80 degrees from flat
- Drop section adjusts from either side and lowers up to 80 degrees from flat
- Smooth, electric, raise and lower height adjustment
- Full width backrest hinge adds rigidity for greater support
- Full width backrest hinge adds rigidity for greater support
- Heavy duty, steel frame is fabricated in the USA by Clinton Industries
- Dependable scissors-lift features an H-leg design for extra stability
- Rounded corner 1" plywood top with 2" of ultra firm padding covered by heavy weight upholstery
- Harden steel bolts and self lubricating pivot bearings for years of dependable service
- Low height may help achieve ADA tax credit


CLINTON OPEN BASE & SHROUDED
POWERXL
Extra width, Extra capacity, Extra Affordable



8039987 SPECIFICATIONS

Overall	76" L x 30" W x 18-35" H
Backrest	38" L x 30" W
Footrest	13" L x 30" W

California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1



Upholstery Colors



OPTIONS

Positioning Straps Set

- 2" x 92" strap is secured with quick hook and loop Velcro® fasteners
- 2 straps per set (92" x 2")

080



Safety Rails (pair)

- 3 height positions
- Plunger height adjustment

082



Hand-Held Controller

- Alternative to foot control
- Stores easily under top when not in use

084



030 Paper Dispenser


- Accommodates paper rolls from 18" (45.72cm) to 21" (53.34cm) wide

030

Paper Cutter

- Clear plastic with elastic strap and buckle

040



Paper Dispenser/Cutter Combo

- Paper dispenser and paper holder ordered as one option

034