

new!

# gen3

surgical case carts



CDS-356-LD shown with optional shelf

### CDS-356-LD-B

Overall outside dimensions:  
29½"W x 28⅛"D x 39"H

Usable Interior Space  
20"W x 24¼"D x 27¼"H

*One shelf configuration recommended*



CDS-356-MS shown with optional shelf

### CDS-356-MS-B

Overall outside dimensions:  
29½"W x 28"D x 39"H

Usable Interior Space  
20"W x 24¼"D x 27¼"H

*One shelf configuration recommended*

## case carts

surgical

### features

- Stainless steel construction.
- Integrated body and floor increases overall cart rigidity.
- High temp, impact resistant bumper mounted directly to cart body reducing the overall footprint.
- Cool-Touch polymer handle.
- 6" diameter stainless steel casters.
- Biohazard Indicator Wheel standard.
- CDS-342 and CDS-345 feature two direction lock casters.
- Roll-out shelves with two position latch mechanism safely locks the shelf in or out.
- Cart is designed to quickly drain during cart washing to minimize residual water.



Cool-Touch Handle



Biohazard Indicator Wheel  
**CartStat**



Modular Bumper



Waterfall Wire Shelf



re•volve™  
caster upgrade



CDS-335 shown with optional shelves

**CDS-333-B**

Overall outside dimensions:  
29½"W x 28"D x 56"H

Usable Interior Space  
20"W x 24¼"D x 44¼"H

*Two shelf configuration recommended*

**CDS-335-B**

Overall outside dimensions:  
35"W x 32¾"D x 56"H

Usable Interior Space  
25½"W x 29"D x 44.¼"H

*Two shelf configuration recommended*

**CDS-342-B**

Overall outside dimensions:  
46¾"W x 28"D x 39"H

Usable Interior Space  
37½"W x 24¼"D x 27"H

*One shelf configuration recommended*

**CDS-345-B**

Overall outside dimensions:  
46¾"W x 28"D x 56"H

Usable Interior Space  
37½"W x 24¼"D x 44¼"H

*Two shelf configuration recommended*

**Shelving Configurations**

**CDS-333-SS** ..... Stainless steel solid shelf for CDS-333, CDS-356-LD and CDS-356-MS.

**CDS-333-SSRO** ... Stainless steel roll-out solid shelf for CDS-333, CDS-356-LD and CDS-356-MS.

**CDS-333-WS** ..... Stainless steel wire shelf for CDS-333, CDS-356-LD and CDS-356-MS.

**CDS-333-WSRO** .. Stainless steel roll-out wire shelf for CDS-333, CDS-356-LD and CDS-356-MS.

**CDS-335-SS** ..... Stainless steel solid shelf for CDS-335.

**CDS-335-SSRO** ... Stainless steel roll-out solid shelf for CDS-335.

**CDS-335-WS** ..... Stainless steel wire shelf for CDS-335.

**CDS-335-WSRO** .. Stainless steel roll-out wire shelf for CDS-335.

**CDS-342-SS** ..... Stainless steel solid shelf for CDS-342 and CDS-345.

**CDS-342-SSRO** ... Stainless steel roll-out solid shelf for CDS-342 and CDS-345.

**CDS-342-WS** ..... Stainless steel wire shelf for CDS-342 and CDS-345.

**CDS-342-WSRO** .. Stainless steel roll-out wire shelf for CDS-342 and CDS-345.

Custom CDS carts engineered to customer specifications available upon request.