

Product Specifications & Standards

- Industry best 3 year warranty
- Weight Capacity: 700 lbs.
- Seat Height: 18" - 37"
ADA Access Board Compliant at the seated and table position
- Length: 60"; extendable to 78"
- SafeGlide™ leg extension length: 18"
- Upholstery width: 28"
- Front storage drawer: 19³/₄" x 19¹/₄" x 4"
- Paper holder accepts 18"-21" rolls
- Optional Mayo Pass-Through Work Surface
- Optional ComfortGrip™ Patient Assist Handles
- 0.9 cubic ft. of storage

Optional Accessories

- 3500687PLUS: EasyGlide™ Mobility System
- 99216: Stainless Steel Treatment Pan
- 99502: One Universal Socket Clamp
- 99504: Deluxe Articulating Knee Crutches*
- 100253: Urology Drain Pan System
- 100407PLUS: Front Drawer Dividers
- 3500077: Footrest Pad
- 100540: Two Support Adapters for Knee Crutches (factory installed)
- 100550: Two Support Adapters for Knee Crutches (field installed)
- 102031PLUS: Brewer Welch Allyn® IV Bracket (factory installed)
- 102032PLUS: Brewer Welch Allyn® IV Bracket (field installed)

* Needs two 99502 and one of either 100540 or 100550

Optional Kits

- 2102894PLUS: Stirrups for 5700PLUS and 5800PLUS (set of two, factory installed)
- 2102898PLUS: Stirrups for 5700PLUS and 5800PLUS (set of two, field installed)
- 2102891PLUS: ComfortGrip™ Patient Assist Handles (set of two, factory installed)
- 2102895PLUS: ComfortGrip™ Patient Assist Handles (set of two, field installed)
- 2102892PLUS: Mayo Pass-Through Work Surface (factory installed)
- 2102896PLUS: Mayo Pass-Through Work Surface (field installed)




EasyGlide™ Mobility System



ComfortGrip™ Patient Assist Handles



Mayo Pass-Through Work Surface



Ergonomic Stirrups

FLEX™ HIGH-LOW PLUS

POWER EXAM TABLE



Optional Features

- EasyGlide™ Mobility System
- ComfortGrip™ Patient Assist Handles
- Ergonomic Stirrups
- Mayo Pass-Through Work Surface

Upholstery Options

- Available in 15 standard colors
- Specialty colors available for an additional cost

Certifications

- UL 60601-1 compliant
- CAN-CSA C22.2 No. 601.1

Additional Specifications for Models 5701PLUS and 5801PLUS

- Hospital-grade outlet - standard
- 115VAC, 60Hz
- Pelvic tilt
- 60 Watt front drawer heater



FLEX High-Low PLUS shown in Blue Fog

Models	Weight Capacity	ADA Compliant 18" Seat Height	Power Base	Power Back	Manual Back	Assist Handles	Stirrups	Pelvic Tilt & Drawer Warmer	Outlet	Work Station	Mobility
5700PLUS	700	•	•		•	optional	optional			optional	optional
5701PLUS	700	•	•		•	optional	•	•	•	optional	optional
5800PLUS	700	•	•	•		optional	optional			optional	optional
5801PLUS	700	•	•	•		optional	•	•	•	optional	optional