

<b>Product Name:</b>	<b>Cover Glasses</b>
<b>Product Application:</b>	<b>Suitable for a variety of cell staining, including immunohistochemistry and immunofluorescence staining followed by laser confocal observation.</b>
<b>Product Materials:</b>	<b>High quality borosilicate glass</b>
<b>product Specifications:</b>	Temperature range: stored at room temperature Shelf life: 3 Years after date of production ( ensure package is in good )
<b>Sterilization:</b>	<b>Yes.</b>
<b>Pyrogens:</b>	Non-pyrogenic
<b>RNase/DNase:</b>	DNase/RNase free
<b>BSE/TSE:</b>	These products are deemed animal free.
<b>Performance Testing:</b>	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: qualified                      Cell culture test: qualified Packaging inspection: qualified
<b>Features</b>	<ol style="list-style-type: none"> <li>1. High quality borosilicate glass, resistant to organic solvents.</li> <li>2. No pyrogen, no cytotoxicity, DNase/RNase free</li> <li>3. Autoclavable</li> <li>4. Not TC treated.</li> <li>5. Thickness: 0.16-0.19mm</li> <li>6. Can be used in confocal experiment.</li> <li>7. Each package has an independent batch number for easy quality tracking and traceability.</li> </ol>
<b>Packages</b>	Each package has a spare quantity to balance the loss during transportation and use.

## Product Range

Cat. No.	Glass Diameter(mm)	Thickness (mm)	Application	/Pack
801010	14	0.16±0.02	24 Well Cell Culture Plate	100
801007	15	0.16±0.02	24 Well Cell Culture Plate	100
801011	18	0.16±0.02	12 Well Cell Culture Plate	100
801008	20	0.16±0.02	12 Well Cell Culture Plate	100
801009	25	0.16±0.02	6 Well Cell Culture Plate	100

Product Picture

